

Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs
and Drug Control

MUSCAT



سِلاطِنَا جَمَانَا
وَزَارَةُ الصِّحَّةِ
وَالرِّقَابَةِ الدَّوَلِيَّةِ
مَسْقَط

To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No. ¹⁷..... dated ^{23/01/20}..... Regarding NCMDR Recall of
Idylla ctNRAS-BRAF Mutation Test from (Mfr: Biocartis)

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information

Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs
and Drug Control
MUSCAT



سِلاطِنَة عُمان
وَزارة الصِّحة
المَدِينَة العَمانَة للصِّدائِكَة
وَالرِّقَابَة الدَّوائِيَّة
مَسَقَط

Circular No. 17 / 2020

Ref: 11 / 2020

27-05-1441 H

23-01-2020

Recall of Idylla ctNRAS-BRAF Mutation Test from Biocartis

| | |
|-----------------------|---|
| Source of Recall | NCMDR National Centre for Medical Devices Reporting, SFDA https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=6&rid=15019 |
| Product | Idylla ctNRAS-BRAF Mutation Test – IVD |
| Manufacturer | Biocartis |
| The affected products | LOT: 00003867, 00003939, 00004216, 00004270 |
| Reason for Recall | After receiving a customer complaint for incorrect labeling, Biocartis has identified that four (4) lots of IdyllarM ctNRAS-BRAF Mutation Test cartridges have been labeled with an incorrect expiry date. |
| Action | <ol style="list-style-type: none">1. Kindly check your stock, contact your local agent for remedial action2. Please stop using all your remaining Idylla™ ctNRAS-BRAF Mutation Test cartridges from the lots listed in Table 1 and destroy them..3. Results generated with impacted Idylla™ ctNRAS-BRAF Mutation Test cartridges, beyond the nine (9) months shelf life and before receipt of this Field Safety Notice, can be considered valid. There is no need to re-test patient samples due to this labelling error. |
| Product image | |
| comments | Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical Device to Director of Medical Device Control contact E-mail dg-padc@moh.gov.om |

Directorate General of Pharmaceutical Affairs & Drug Control
Sultanate of Oman

Dr. Mohammed Hamdan Al Rubaie

DIRECTOR GENERAL



For the Attention of the Laboratory Director

URGENT – Field Safety Notice

Idylla™ ctNRAS-BRAF Mutation Test Issue: Incorrect expiry date on product label

| | |
|---|--|
| Product Name | Idylla™ ctNRAS-BRAF Mutation Test |
| Device Identifier | |
| REF | A0090/6 |
| GTIN | 15415219000710 |
| Production Identifier (Lot. No.) | 00003867 00003939 00004216 00004270 |
| Type of Action | Field Safety Corrective Action |

Dear Valued Customer,

Biocartis has identified that the labeling of the above listed lots (see Production Identifier) of Idylla™ ctNRAS-BRAF Mutation Test cartridges contains an incorrect expiry date. Due to this labeling error, Biocartis has initiated a Field Safety Corrective Action to prevent further use of the impacted Idylla™ ctNRAS-BRAF Mutation Test cartridges.

Problem Description

After receiving a customer complaint for incorrect labeling, Biocartis has identified that four (4) lots of Idylla™ ctNRAS-BRAF Mutation Test cartridges have been labeled with an incorrect expiry date. The shelf life of the concerned product was erroneously set to twelve (12) months instead of to the currently claimed shelf life of nine (9) months. The impacted lots are listed in Table 1 below, together with the incorrect expiry dates (as marked on the product label) and the actual (correct) expiry dates.

Table 1: Overview expiry dates of impacted Idylla™ ctNRAS-BRAF Mutation Test cartridges

| Lot. No. | Expiry date on label (incorrect) | Correct expiry date |
|----------|----------------------------------|---------------------|
| 00003867 | 2019-07-18 | 2019-04-14 |
| 00003939 | 2019-12-06 | 2019-09-02 |
| 00004216 | 2020-05-29 | 2020-02-24 |
| 00004270 | 2020-07-08 | 2020-04-04 |

Potential risk

Two of the affected lots (00003867, 00003939), which were distributed to the market, could potentially have been used beyond their claimed shelf life of nine (9) months (see Table 1).

Although the claimed shelf life of the Idylla™ ctNRAS-BRAF Mutation Test has not been extended to twelve (12) months, available stability data demonstrates that the product performance remains within acceptance criteria after twelve (12) months of storage. In addition, Biocartis has tested performance of the retained cartridges from the oldest impacted lot (00003867) after thirteen (13) months of storage. All tested cartridges provided a correct result (i.e. correct call). This indicates that there is no increased risk of generating incorrect results with the impacted Idylla™ ctNRAS-BRAF Mutation Test cartridges, when used between nine (9) and twelve (12) months of storage.

Actions to be taken by the customer

- 1) Please stop using **all** your remaining Idylla™ ctNRAS-BRAF Mutation Test cartridges from the lots listed in Table 1 and destroy them.
- 2) Please return the completed 'Acknowledgement of Receipt' form in Appendix 1 of this Field Safety Notice to Biocartis to confirm the destruction of the impacted Idylla™ ctNRAS-BRAF Mutation Test cartridges.
- 3) Results generated with impacted Idylla™ ctNRAS-BRAF Mutation Test cartridges, beyond the nine (9) months shelf life and before receipt of this Field Safety Notice, can be considered valid. There is no need to re-test patient samples due to this labeling error.

Actions taken by Biocartis NV

- 1) Biocartis has notified local Regulatory Authorities of this Field Safety Corrective Action.
- 2) The root cause of the event has been identified, and immediate actions have been taken to prevent that any remaining non-conforming product is distributed to the market. Implementation of corrective actions is ongoing to prevent recurrence of the event.
- 3) Biocartis will replace, free of charge, all Idylla™ ctNRAS-BRAF Mutation Test cartridges for which Biocartis receives confirmation of destruction through means of the completed Appendix 1 of this Field Safety Notice.

Biocartis Field Safety Notice
Biocartis Reference: BC-012798 Rev. 1
Date: September 11, 2019



Communication of this Field Safety Notice

Please forward this information to all individuals and departments within your organization that have received or used this product. If you are not the end user, please forward this Field Safety Notice to the device end user. Please maintain the awareness of this Field Safety Notice and resulting action for an appropriate period to ensure the effectiveness of the corrective action.

Completion of the 'Acknowledgement of Receipt' form

Due to regulatory reasons, completion of the 'Acknowledgement of Receipt' form (Appendix 1 of this Field Safety Notice) is required. Please complete and sign the attached 'Acknowledgement of Receipt' form by September 20, 2019 and email it to hotline@biocartis.com.

We sincerely apologize for any inconvenience this may cause and thank you in advance for your understanding and support.

If you need any further information or assistance concerning this notice, please contact the Biocartis hotline (Phone: +32 (0) 15 632 800 between 9h00 and 17h00 CEST; e-mail: hotline@biocartis.com) or your local Biocartis representative.

Yours sincerely,



Biocartis NV