



نتقدم بثقة  
Moving Forward  
with Confidence

رؤية عُمان  
2040  
Oman Vision

To:

**THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES**

**Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)**

**Director General of Engineering Affairs, MOH**

**Director General of Royal Hospital**

**Director General of Khoula Hospital**

**Director General of Medical Supplies (MOH)**

**Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)**

**Hospital Director (Al Nahda Hospital)**

**Hospital Director (Al Massara Hospital)**

**The Head of Medical Services in SQU Hospital**

**The Head of Medical Services in Royal Oman Police**

**The Head of Medical Services in Ministry of Defence**

**The Head of Medical Services in The Diwan**

**The Head of Medical Services in The Sultan's Special Force**

**The Head of Medical Services in Internal Security Services**

**The Head of Medical Services in Petroleum Development of Oman**

**The Head of Medical Services in LNG Oman**

**ALL PRIVATE PHARMACIES & DRUG STORES**

After Compliments,

Please find attached our Circular No **66** dated **19/5/2024** Regarding NCMDR Field Safety Corrective Action of Various Mapleson F Anaesthetic Breathing Systems from (mfr: Intersurgical Limited).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DSC
- Director of Pharmacovigilance & Drug Information Dept, DSC
- Director of Drug Control Department, DSC
- Director of Pharmaceutical Licensing Department, DSC
- Director of Central Quality Control Lab., DSC
- Supdt. of Central Drug Information



**DSC**  
مركز سلامة الدواء  
Drug Safety Center



ص.ب: ٣٩٣ مسقط - الرمز البريدي: ١٠٠ - هاتف: ٢٢٣٥٧١١١ - فاكس: ٢٢٣٥٨٤٨٩

P.O. Box: 393 Muscat - Postal Code: 100 - Tel: 22357111 - Fax: 22358489

✉ @DSCPHO Email: dscpho@moh.gov.om



Circular No. 66 / 2024

11-11-1445 H  
19-05-2024

نتقدم بثقة  
Moving Forward  
with Confidence

رؤية عمان  
2040  
Oman Vision

Field Safety Corrective Action of Various Mapleson F Anaesthetic Breathing Systems from Intersurgical Limited.

Source	NCMDR - National Center Medical Device Reporting- SFDA. <a href="https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=21038">https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=21038</a>
Product	Various Mapleson F Anaesthetic Breathing Systems.
Description	Breathing circuit, anaesthesia, single use.
Manufacturer	Intersurgical Limited.
Local agent	Muscat Pharmacy & Stores LLC.
The affected products	Device Model/Catalogue/part number(s): 2120000, 2121000, 2121002, 2121004, 2121005, 2121011, 2121014, 2121019, 2121024, 2121035, 2121042, 2121045, 2121048, 2121053 Expiry date from April 2024 to March 2029 Please refer to the attachment for more information.
Reason	Some devices contain a reservoir bag with a closed tail, when they should have a reservoir bag with open tail.
Action	1. Identify any potentially affected products from the affected codes and lot numbers listed in the attachment. 2. All users must perform a thorough visual inspection and functional test before use of the affected products, to confirm a patent gas pathway exists through the open tail of the reservoir bag to avoid over pressurisation of the system. 3. Retain any affected sample(s) identified. 4. The manufacturer will also be introducing a new instruction for use which will include the following Pre-Use Check in line with the recommended action mentioned above: 5. "If the product is supplied without an APL valve, the pressure within the system is controlled by the clinician through manipulation of the open tail on the reservoir bag. Check that a patent gas pathway exists through the open tail of the reservoir bag". 6. Contact the local agent for remedial action.
comments	Healthcare professionals are encouraged to report any adverse events suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: <a href="mailto:Med-device@moh.gov.om">Med-device@moh.gov.om</a>

Dr. Mohammed Hamdan Al Rubaie  
Director General



Date: 24.04.2024

## **Urgent Field Safety Notice**

### **VARIOUS MAPLESON F ANAESTHETIC BREATHING SYSTEMS CONTAINING 0.5L RESERVOIR BAGS WITH OPEN TAILS**

For Attention of\*: MDSO's, All clinical staff, Managers and users of the above product

Contact details of local representative (name, e-mail, telephone, address etc.)\*

**Giedrius Budrys**  
**Customer Resolution and Relationship Manager**  
**Intersurgical UAB**  
**Arnioniu str 60, LT-18170 Pabrade Lithuania**

**Email: [giedriusb@intersurgical.lt](mailto:giedriusb@intersurgical.lt)**  
**Tel. +370 387 66611**  
**Fax: +370 387 66622**

or

**This could be a distributor or local branch of the manufacturer. To be added at  
the appropriate stage in the different local languages**





## **Urgent Field Safety Notice (FSN)**

### **VARIOUS MAPLESON F ANAESTHETIC BREATHING SYSTEMS CONTAINING 0.5L RESERVOIR BAGS WITH OPEN TAILS**

#### **Risk addressed by FSN**



<b>1. Information on Affected Devices*</b>																			
1	1. Device Type(s)*																		
.	Various Mapleson F Anaesthetic Breathing Systems																		
1	2. Commercial name(s)																		
.	Mapleson F infant T-piece breathing system with 0.5L open tail bag, $\geq 1.8\text{m}$ Mapleson F Jackson Rees modification T-piece breathing system with 0.5L open tail bag, $\geq 1.8\text{m}$ Paediatric, Mapleson F Jackson Rees modification T-piece with 0.5L open tail bag, and swivel elbow, $\geq 1.8\text{m}$ Paediatric, Mapleson F Jackson Rees modification T-piece with 0.5L open tail bag, and swivel elbow, $\geq 4.8\text{m}$ Paediatric, Mapleson F Jackson Rees modification T-piece with 0.5L open tail bag, and elbow, $\geq 2.8\text{m}$ Paediatric, Mapleson F Jackson Rees modification T-piece with 0.5L open tail bag, and elbow, $\geq 1.8\text{m}$ <b>Paediatric, Mapleson F Jackson Rees modification T-piece with 0.5L open tail bag, and swivel elbow, <math>\geq 3.6\text{m}</math></b> Paediatric, Mapleson F Jackson Rees modification T-piece with 0.5L open tail bag, and luer elbow, $\geq 3.6\text{m}$ Map/F 0.5L Open T/B Luer/Elb $\geq 2.4\text{m}$ Map/F 0.5L Open T/B Luer/Elb M/Line $\geq 1.8\text{m}$ Map/F 0.5L Open T/B Luer/Elb $\geq 1.6\text{m}$ Paediatric, Mapleson F Jackson Rees modification T-piece with 0.5L open tail bag, and luer elbow, $\geq 10.8\text{m}$																		
1	3. Unique Device Identifier(s) (UDI-DI)																		
.	<table border="0" style="width: 100%;"> <tr> <td>5030267062249</td><td>5030267062270</td><td>5030267062362</td><td>5030267062430</td><td>5030267103164</td></tr> <tr> <td>5030267062256</td><td>5030267062287</td><td>5030267062379</td><td>5030267062508</td><td>5030267149810</td></tr> <tr> <td>5030267062263</td><td>5030267062348</td><td>5030267062393</td><td>5030267062539</td><td></td></tr> </table>				5030267062249	5030267062270	5030267062362	5030267062430	5030267103164	5030267062256	5030267062287	5030267062379	5030267062508	5030267149810	5030267062263	5030267062348	5030267062393	5030267062539	
5030267062249	5030267062270	5030267062362	5030267062430	5030267103164															
5030267062256	5030267062287	5030267062379	5030267062508	5030267149810															
5030267062263	5030267062348	5030267062393	5030267062539																
	4. Primary clinical purpose of device(s)*																		
	To deliver and remove anaesthetic and respiratory gases to and from a paediatric patient via a breathing system comprised of tubing and connectors and 0.5 L reservoir bag.																		
1	5. Device Model/Catalogue/part number(s)*																		
.	2120000, 2121000, 2121002, 2121004, 2121005, 2121011, 2121014, 2121019, 2121024, 2121035, 2121042, 2121045, 2121048, 2121053																		
1	6. Software version																		
.	N/A																		
1	7. Affected serial or lot number range																		
.	Any of the above with an expiry date from April 2024 to March 2029.																		
1	8. Associated devices																		
.	N/A.																		

<b>2. Reason for Field Safety Corrective Action (FSCA)*</b>	
2.	<p>1. Description of the product problem*</p> <p>Some devices contain a reservoir bag with a closed tail, when they should have a reservoir bag with open tail.</p> <div style="display: flex; justify-content: space-around; align-items: center;">   </div> <div style="display: flex; justify-content: space-around; align-items: center;"> <p><b>Correct – Open Tail Reservoir Bag</b></p> <p><b>Incorrect – Closed Tail Reservoir Bag</b></p> </div>
2.	<p>2. Hazard giving rise to the FSCA*</p> <p>If the incorrect closed tail reservoir bag is not identified during the routine pre use check as described in the product instruction for use, it could result in over pressurisation of the system leading to potential barotrauma.</p>
2.	<p>3. Probability of problem arising</p> <p>100% in the affected range.</p>
2.	<p>4. Predicted risk to patient/users</p> <p>The risks associated with the identified fault have been reviewed, and whilst the probability of occurrence is low, we believe it is essential to address the issue promptly to further reduce the risk of any potential patient harm.</p>
2.	<p>5. Further information to help characterise the problem</p> <p>N/A</p>
2.	<p>6. Background on Issue</p> <p>Following a customer report from the market and subsequent thorough inspection and analysis of internal stock, we have identified a potential safety concern related to various Mapleson F paediatric anaesthetic breathing systems as listed above. Unfortunately some products have been manufactured with a 0.5 L reservoir bag with a closed tail which could result in over pressurisation of the system.</p>
2.	<p>7. Other information relevant to FSCA</p> <p>N/A</p>
	<b>3. Type of Action to mitigate the risk*</b>



<b>3.</b>	<b>1. Action To Be Taken by the User*</b>	
	<div style="display: flex; justify-content: space-between;"> <span><input checked="" type="checkbox"/> Identify Device</span> <span><input type="checkbox"/> Quarantine Device</span> <span><input type="checkbox"/> Return Device</span> <span><input type="checkbox"/> Destroy Device</span> </div> <div style="display: flex; justify-content: space-between;"> <span><input type="checkbox"/> On-site device modification/inspection</span> </div> <div style="display: flex; justify-content: space-between;"> <span><input type="checkbox"/> Follow patient management recommendations</span> </div> <div style="display: flex; justify-content: space-between;"> <span><input checked="" type="checkbox"/> Take note of amendment/reinforcement of Instructions For Use (IFU)</span> </div> <div style="display: flex; justify-content: space-between;"> <span><input checked="" type="checkbox"/> Other</span> <span><input type="checkbox"/> None</span> </div> <p>Please distribute this Field Safety Notice to all potential users of the Mapleson F paediatric anaesthetic breathing systems listed above, within your facility. This is for their awareness of the potential problem and to carry out the following actions.</p> <p>To ensure the safety of patients we recommend the following actions.</p> <ol style="list-style-type: none"> <li>1. Identify any potentially affected products from the affected codes and lot numbers listed above.</li> <li>2. All users must perform a thorough visual inspection and functional test before use of the products and lot numbers listed above, to confirm a patent gas pathway exists through the open tail of the reservoir bag to avoid over pressurisation of the system.</li> <li>3. Retain any affected sample(s) identified, and please report to us immediately.</li> </ol> <p><b>Please note:</b> This is not a product removal.</p> <p>Please complete and return the Reply Form provided to <a href="mailto:giedriusb@intersurgical.lt">giedriusb@intersurgical.lt</a> (local contact e-mail address), to confirm receipt of this notice and that the necessary actions are being taken.</p> <p>Please continue to report to Intersurgical any adverse events involving this product.</p>	
<b>3.</b>	<b>2. By when should the action be completed?</b>	Immediately on receipt of this FSN, and awareness of this FSN should be ongoing until all potentially affected stock listed in this FSN has been used up.
<b>3.</b>	<b>3. Particular considerations for: N/A</b>  Is follow-up of patients or review of patients' previous results recommended?  Not applicable.	
<b>3.</b>	<b>4. Is customer Reply Required? *</b> (If yes, form attached specifying deadline for return)	Yes
<b>3.</b>	<b>5. Action Being Taken by the Manufacturer</b>  <div style="display: flex; justify-content: space-between;"> <span><input type="checkbox"/> Product Removal</span> <span><input type="checkbox"/> On-site device modification/inspection</span> </div> <div style="display: flex; justify-content: space-between;"> <span><input type="checkbox"/> Software upgrade</span> <span><input checked="" type="checkbox"/> IFU or labelling change</span> </div> <div style="display: flex; justify-content: space-between;"> <span><input checked="" type="checkbox"/> Other</span> <span><input type="checkbox"/> None</span> </div>	

	<p>We have implemented corrective actions in manufacturing process to eliminate this problem for future supply. We will also be introducing a new instruction for use which will include the following Pre-Use Check in line with the recommended action 2. above:</p> <p><i>If the product is supplied without an APL valve, the pressure within the system is controlled by the clinician through manipulation of the open tail on the reservoir bag. Check that a patent gas pathway exists through the open tail of the reservoir bag.</i></p>	
3	6. By when should the action be completed?	One month from receipt of the FSN
3.	7. Is the FSN required to be communicated to the patient /lay user?	No
3	8. If yes, has manufacturer provided additional information suitable for the patient/lay user in a patient/lay or non-professional user information letter/sheet?	
	N/A	

<b>4. General Information*</b>		
4.	1. FSN Type*	New – Advisory Notice
4.	2. For updated FSN, reference number and date of previous FSN	N/A
4.	3. For Updated FSN, key new information as follows:	
	N/A	
4.	4. Further advice or information already expected in follow-up FSN? *	No
4	5. If follow-up FSN expected, what is the further advice expected to relate to:	
	N/A	
4	6. Anticipated timescale for follow-up FSN	N/A
4.	7. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	Intersurgical Ltd.
	b. Address	Crane House, Molly Millars Lane, Wokingham, Berkshire, RG41 2RZ
	c. Website address	<a href="https://www.intersurgical.com/">https://www.intersurgical.com/</a>
4.	8. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *	
4.	9. List of attachments/appendices:	Customer Reply Form
4.	10. Name/Signature	Ivan Seniut, Group Quality and Regulatory Affairs Director, Intersurgical
	<p>E-Signed by Ivan Seniut VERIFY authenticity with ApproveIt </p> 	

Transmission of this Field Safety Notice	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to <b>ensure effectiveness of the corrective action.</b></p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.</p>

Note: Fields indicated by \* are considered necessary for all FSNs. Others are optional.