



To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES
Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)
Director General of Engineering Affairs, MOH
Director General of Royal Hospital
Director General of Khoula Hospital
Director General of Medical Supplies (MOH)
Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)
Hospital Director (Al Nahda Hospital)
Hospital Director (Al Massara Hospital)
The Head of Medical Services in SQU Hospital
The Head of Medical Services in Royal Oman Police
The Head of Medical Services in Ministry of Defence
The Head of Medical Services in The Diwan
The Head of Medical Services in The Sultan's Special Force
The Head of Medical Services in Internal Security Services
The Head of Medical Services in Petroleum Development of Oman
The Head of Medical Services in LNG Oman
ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No. 279 dated 25/12/23 Regarding NCMDR Field Safety Corrective Action of Neuroinspire neurolocate - Neuromate Stereotactic System from (mfr: Renishaw Mayfield).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information





Circular No. 279/2023

11 -06-1445 H

25 -12-2023

لتقدم بثقة
Moving Forward
With Confidence

رؤية عمان
2040
Oman Vision

**Field Safety Corrective action of Neuroinspire neurolocate - Neuromate Stereotactic System
from Renishaw Mayfield**

Source	NCMDR- National Center for Medical Devices Reporting- SFDA https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&rid=19830
Product	Neuroinspire neurolocate - Neuromate Stereotactic System.
Description	Neurosurgical robot for stereotactic applications, used for the positioning of a tool guide.
Manufacturer	Renishaw Mayfield.
The affected products	Software version: V6.2.2 and above IFU Version number is found on the front page as shown in the examples in the attachment. Model/Catalogue/part number of the IFU: 047.0107C/EN neuroinspire neurolocate module – English Language.
Reason	Some of the images in the IFUs are incorrect and do not relate to the corresponding text in the IFUs.
Action	1. The images in the IFUs cited in section 1 in the attachment will be updated and new IFUs provided to affected sites. 2. Contact the local agent for remedial action.
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: Med-device@moh.gov.om

Dr. Mohammed Hamdan Al Rubaie
Director General



PADC
المديرية العامة للصيدلة والرقابة الدوائية
Directorate General of Pharmaceutical
Affairs & Drug Control



ص.ب: 393 مسقط - الرمز البريدي: 100 - هاتف: 22357111 - فاكس: 22358489
P.O. Box: 393 Muscat - Postal Code: 100 - Tel: 22357111 - Fax: 22358489
dgpa_dc Email: dg-padc@moh.gov.om

FSN Ref: VR-23-12.1

FSCA Ref: VR-23-12

Date: 2023-12-08

Urgent Field Safety Notice
Neuromate Stereotactic System
Neuroinspire – neurolocate module
Incorrect images in the IFU





For Attention of: Neurosurgeons, Neurosurgical Theatre Staff, Hospital Device Safety Officer

Contact details of local representative (name, e-mail, telephone, address etc.)*

R Rusling
Renishaw Mayfield SARL
31 rue Ampère
Chassieu
69680
France

RNSRegulatory@renishaw.com

Urgent Field Safety Notice
Neuromate Stereotactic System
Neuroinspire – neurolocate module
Incorrect images in the IFU

1. Information on Affected Devices*	
1.	<p>1. Device Type(s)*</p> <p>IFU for neuroinspire neurolocate module within neuroinspire surgical planning software (the software for the neuromate stereotactic robot)</p>
1.	<p>2. Commercial name(s)*</p> <p>neuroinspire neurolocate (optional module in the neuroinspire software)</p>
1.	<p>3. Unique Device Identifier(s) (UDI-DI)</p> <p>Not Applicable for IFU</p>
1.	<p>4. Primary clinical purpose of device(s)*</p> <p>The device is used for stereotactic brain surgery. The affected part is the IFU for the additional/optional module for neuroinspire neurolocate .</p>
1.	<p>5. Device Model/Catalogue/part number(s)*</p> <p>047.0107C/DE neuroinspire neurolocate module - German Language 047.0107B/ES neuroinspire neurolocate module – Spanish Language 047.0107C/FR neuroinspire neurolocate module – French Language 047.0107C/IT neuroinspire neurolocate module – Italian Language 047.0107C/EN neuroinspire neurolocate module – English Language</p> <p>IFU Version number is found on the front page as shown in the examples below:</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 2px; font-size: 8px;"> neuroinspire™ neurolocate™-Modul – Bedienungsanleitung 047.0107C/DE </div> <div style="text-align: center;">  </div> </div> <div style="margin: 10px 0;">  <p style="text-align: center; font-size: 8px;">zur Verwendung mit dem neurolocate™-Modul</p> </div> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="border: 1px solid black; padding: 2px; font-size: 8px;"> neuroinspire™ neurolocate™ module – instructions for use 047.0107C/EN </div> <div style="text-align: center;">  </div> </div> <div style="margin: 10px 0;">  <p style="text-align: center; font-size: 8px;">for use with neurolocate™ module</p> </div>

FSN Ref: VR-23-12.1

FSCA Ref: VR-23-12

1.	6. Software version
	V6.2.2 and above
1.	7. Affected serial or lot number range
	IFUs as detailed above
1.	8. Associated devices
	N/A

2. Reason for Field Safety Corrective Action (FSCA)*	
2.	<p>1. Description of the product problem*</p> <p>Some of the images in the IFUs are incorrect and do not relate to the corresponding text in the IFUs.</p> <p><u>Error 1</u></p> <p>047.0107C/DE neuroinspire neurolocate module - German Language 047.0107B/ES neuroinspire neurolocate module – Spanish Language 047.0107C/FR neuroinspire neurolocate module – French Language 047.0107C/IT neuroinspire neurolocate module – Italian Language</p> <p>The images are consistently incorrect across all languages noted above.</p> <p>Section 3.8 – Figure 27 Section 3.9 – Figure 28</p> <p>There is no defect with the neuroinspire software or the neuromate robot.</p> <hr/> <p>Some of the images in the IFUs are incorrect and do not relate to the corresponding text in the IFUs.</p> <p><u>Error 2</u></p> <p>047.0107C/EN – neuroinspire neurolocate module – English Language</p> <p>Page 1 – Front Image Section 2.2 – Figure 1 Section 2.4 – Figure 3 Section 2.5 – Figure 4 Section 2.6 – Figure 5 Section 2.7 – Figure 6 Section 2.7 – Figure 7 Section 3.3 – Figure 22 Section 3.4 – Figure 23 Section 3.5 – Figure 24 Section 3.6 – Figure 25 Section 3.7 – Figure 26 Section 3.8 – Figure 27 Section 3.8 – Figure 28</p>

	<p>Unless noted in error 1, the above images in error 2 only impact the English language version of the IFU.</p> <p>There is no defect with the neuroinspire software or the neuromate robot.</p>
2.	<p>2. Hazard giving rise to the FSCA*</p> <p>There is likely to be confusion if the surgeon or theatre staff try to use the images to support the text in the IFUs stated above.</p>
2.	<p>3. Probability of problem arising</p> <p>The neurolocate workflow has been designed to be very prescriptive using a step by step process within the software, the incorrect images presented in the IFUs are not available to select within the workflow as an option. These image errors would be very obvious to a trained user, they would be expected to be able to identify these errors. Therefore, a problem arising from these image errors is extremely unlikely.</p> <p>The images are obviously wrong and make no sense in relation to the corresponding text. The text is correct and is sufficient to support the trained user in the use of the software. All users are fully trained in the use of the system by Renishaw Field Service personnel.</p> <p>Error 1</p> <p>Figures 27 & 28 are images which are neurolocate module specific. The text is clear on the actions required, the images should give the expected view on the screen. The software user interface and the text in the IFU are sufficient for the surgeon to continue the surgery with no harm to the patient. The correct images are presented below and provide the correct view on the screen</p> <p>Error 2</p> <p>All the Figures are images which are neurolocate module specific. The text is clear on the actions required, the images give the expected view on the screen. The software user interface and the text in the IFU are sufficient for the surgeon to continue the surgery with no harm to the patient.</p>
2.	<p>4. Predicted risk to patient/users</p> <p>Very unlikely to result in patient harm</p>
2.	<p>5. Further information to help characterise the problem</p> <p>None</p>
2.	<p>6. Background on Issue</p> <p>This issue was detected within the Renishaw organisation. It has not been identified by a clinician in the field.</p>
2.	<p>7. Other information relevant to FSCA</p>
	<p>Previous Versions of the IFUs Impacted</p>

The errors presented also impacted historic versions of this IFU. The identification of 1 of the historic IFUs follow a different numeric reference type, this was then subsequently transferred to a new numeric reference type.

Error 1

- H-5630-3035-A3 – neurolocate module – instructions for use (English Language, French Language, Spanish Language, German Language, Italian Language)
- 047.0107A/EN – neurolocate module – instructions for use (English Language)
- 047.0107B/DE – neurolocate module – instruction for use (German Language)
- 047.0107B/FR – neurolocate module – instructions for use (French Language)
- 047.0107B/IT – neurolocate module – instructions for us (Italian Language)





Please note that version A was not created for German, French, Italian or Spanish language.

Error 2

- H-5630-3035-A3 – neuroinspire surgical planning software for use with neurolocate module
- 047.0107A/EN – neurolocate module – instructions for use (English Language)
- 047.0107B/EN – neurolocate module – instructions for use (English Language)







Details of the correct and incorrect images are detailed below :

Error 1

Correct ✓	Incorrect ✗
<p>Figure 27</p> 	<p>Figure 27</p> 
<p>Figure 28</p> 	<p>Figure 28</p> 

Error 2

The incorrect images where there is most likely to be confusion are detailed below.

Correct ✓	Incorrect ✗
<p data-bbox="496 554 823 646">2.5 neurolocate 2D patient localization process Figure 4</p> 	<p data-bbox="1038 554 1366 646">2.5 neurolocate 2D patient localization process Figure 4</p> 
<p data-bbox="560 919 679 952">Figure 27</p> 	<p data-bbox="1110 919 1230 952">Figure 27</p> 
<p data-bbox="560 1247 679 1279">Figure 28</p> 	<p data-bbox="1110 1247 1230 1279">Figure 28</p> 

A full list of images for Error 2 can be supplied upon request.

1. Type of Action to mitigate the risk*	
3.	<p>1. Action To Be Taken by the User*</p> <p><input checked="" type="checkbox"/> Identify Device <input type="checkbox"/> Quarantine Device <input type="checkbox"/> Return Device <input type="checkbox"/> Destroy Device</p> <p><input type="checkbox"/> On-site device modification / inspection</p> <p><input type="checkbox"/> Follow patient management recommendations</p> <p><input type="checkbox"/> Take note of amendment / reinforcement of Instructions For Use (IFU)</p> <p><input checked="" type="checkbox"/> Other <input type="checkbox"/> None</p> <p>Provide further details of the action(s) identified.</p>
3.	<p>2. By when should the action be completed? Review should be immediate to identify incorrect IFUs</p>
3.	<p>3. Particular considerations for: Choose an item.</p> <p>Is follow-up of patients or review of patients' previous results recommended? No</p> <p>No impact on patient expected</p>
3.	<p>4. Is customer Reply Required? * No (If yes, form attached specifying deadline for return)</p>
3.	<p>5. Action Being Taken by the Manufacturer*</p> <p><input type="checkbox"/> Product Removal <input type="checkbox"/> On-site device modification/inspection <input type="checkbox"/> Software upgrade <input checked="" type="checkbox"/> IFU or labelling change <input type="checkbox"/> Other <input type="checkbox"/> None</p> <p>The images in the IFUs cited in section 1 will be updated and new IFUs provided to affected sites.</p>
3.	<p>6. By when should the action be completed? January 2024</p>
3.	<p>7. Is the FSN required to be communicated to the patient /lay user? No</p>
3.	<p>8. If yes, has manufacturer provided additional information suitable for the patient/lay user in a patient/lay or non-professional user information letter/sheet?</p> <p>Choose an item. Choose an item.</p>

2. General Information*		
4.	1. FSN Type*	New
4.	2. For updated FSN, reference number and date of previous FSN	Provide reference and date of previous FSN if relevant.
4.	3. For Updated FSN, key new information as follows: Summarise any key difference in devices affected and/or action to be taken.	
4.	4. Further advice or information already expected in follow-up FSN? *	No
4.	5. If follow-up FSN expected, what is the further advice expected to relate to: Eg patient management, device modifications etc.	
4.	6. Anticipated timescale for follow-up FSN	For provision of updated advice.
4.	7. Manufacturer information (For contact details of local representative refer to page 1 of this FSN)	
	a. Company Name	See Page 1
	b. Address	See Page 1
	c. Website address	See Page 1
4.	8. The Competent (Regulatory) Authority of your country has been informed about this communication to customers. *	
4.	9. List of attachments/appendices:	If extensive consider providing web-link instead.
4.	10. Name/Signature	R Rusling Regulatory Affairs Manager

Transmission of this Field Safety Notice	
	<p>This notice needs to be passed on all those who need to be aware within your organisation or to any organisation where the potentially affected devices have been transferred. (As appropriate)</p> <p>Please transfer this notice to other organisations on which this action has an impact. (As appropriate)</p> <p>Please maintain awareness on this notice and resulting action for an appropriate period to ensure effectiveness of the corrective action.</p> <p>Please report all device-related incidents to the manufacturer, distributor or local representative, and the national Competent Authority if appropriate, as this provides important feedback.*</p>

Note: Fields indicated by * are considered necessary for all FSNs. Others are optional.