



To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES
Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)
Director General of Engineering Affairs, MOH
Director General of Royal Hospital
Director General of Khoula Hospital
Director General of Medical Supplies (MOH)
Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)
Hospital Director (Al Nahda Hospital)
Hospital Director (Al Massara Hospital)
The Head of Medical Services in SQU Hospital
The Head of Medical Services in Royal Oman Police
The Head of Medical Services in Ministry of Defence
The Head of Medical Services in The Diwan
The Head of Medical Services in The Sultan's Special Force
The Head of Medical Services in Internal Security Services
The Head of Medical Services in Petroleum Development of Oman
The Head of Medical Services in LNG Oman
ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No 274 dated 21/12/2023 Regarding NCMDR Field Safety Notice of SynchroMed™ II from (mfr: Medtronic Inc).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information





بتقدم بثقة
Moving Forward
with Confidence



Circular No. 274 / 2023

08 -06-1445 H

21 -12-2023

Field Safety Notice of SynchroMed™ II from Medtronic Inc.

Source	NCMDR - National Center Medical Device Reporting- SFDA. https://ncmdr.sfda.gov.sa/Secure/CA/CAViewRecall.aspx?caid=4&rid=19807
Product	SynchroMed™ II.
Description	Implantable Programmable Infusion Pump.
Manufacturer	Medtronic Inc.
Local agent	Taiba Medserve.
The affected products	Model 8637 Product Number/ CFN: PUMP 8637-20 SYNCHROMED II (8637-20), UDI-Device Identifier (GTIN/UPN): 00643169345201, 00643169384101, 00643169700925, 00643169700932, 00643169732278, 00643169999862, 00763000122676, 00763000421748, 00763000421755, 00763000689643 Product Number/ CFN: PUMP 8637-40 SYNCHROMED II (8637-40), UDI-Device Identifier (GTIN/UPN): 00643169345232, 00643169384170, 00643169384231, 00643169701021, 00643169701038, 00643169701090, 00643169732377, 00643169999961, 00763000122768, 00763000421830, 00763000421847, 00763000689575
Reason	When turning on BIOFIRE TORCH systems, arcing inside of the power switch ultimately cause deformation may result in an open circuit causing the power switch to fail. This event would only occur after the product is in use.
Action	1. It has been recently identified by Medtronic that if the SynchroMed II pump switches into telemetry mode due to electromagnetic interference (EMI) from an MRI scan, while the pump is sounding an alarm, the pump will not resume drug delivery after leaving the MRI magnetic field. 2. Please follow "Patient Management Recommendations" in the attachment. 3. Contact the local agent for remedial action.
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: Med-device@moh.gov.om

Dr. Mohammed Hamdan Al Rubaie

Director General



PADC
المديرية العامة للصيدلة والرقابة الدوائية
Directorate General of Pharmaceutical
Affairs & Drug Control



ص.ب. 393 مسقط - الرمز البريدي 22357111 - فاكس: 22358489
P.O. Box: 393 Muscat - Postal Code: 100 - Tel: 22357111 - Fax: 22358489
dgpa_dc Email: dg-padc@moh.gov.om



Urgent Field Safety Notice

Model 8637 SynchroMed™ II

MRI Guidelines for the SynchroMed Infusion System

November 2023

Medtronic Reference: FA1367

Dear Healthcare Professional,

The purpose of this letter is to communicate the need to interrogate the SynchroMed II™ pump following Magnetic Resonance Imaging (MRI).

Issue Description:

The *MRI Guidelines Instructions for Use for Medtronic Model 8637 Implantable Infusion Systems* (MRI Guidelines) indicate that during normal operations, the magnetic field of the MRI scanner will temporarily stop the rotor of the SynchroMed II pump motor and suspend drug infusion for the duration of the MRI exposure. The pump should resume normal operation upon termination of MRI exposure.

Medtronic recently identified that if the SynchroMed II pump switches into telemetry mode due to electromagnetic interference (EMI) from an MRI scan, while the pump is sounding an alarm, the pump will **not resume** drug delivery after leaving the MRI magnetic field, which is inconsistent with the current labeling. In this case, drug delivery will only resume after performing a post-MRI pump interrogation with the Clinician Programmer (or Personal Therapy Manager) which will end telemetry mode.

If the SynchroMed II pump does not resume drug delivery after leaving the MRI magnetic field, patients may experience a return of underlying symptoms (i.e., pain or spasticity) due to loss of therapy, potentially requiring outpatient or inpatient management, and in severe cases (i.e., baclofen withdrawal), life-threatening or fatal withdrawal symptoms could occur.

From January 01 2019 through October 18 2023, Medtronic has received a total of 13 complaints related to this issue. The complaints reported non-serious underdose symptoms (i.e., withdrawal or return of symptoms) when a follow-up interrogation was not performed post-MRI. After the pump was interrogated, the issue was resolved, and therapy resumed.

Refer to Appendix 1 for product scope.

Medtronic

Patient Management Recommendations:

- Upon completion of an MRI scan, interrogate the pump with the Clinician Programmer (or Personal Therapy Manager) to end telemetry mode and resume drug delivery.
- Consult the MRI Guidelines for additional information on MRI preparation and post-examination review, and motor stall recovery timing (see MRI Guidelines at www.manuals.medtronic.com).
- Remind your patients about the importance of interrogating the SynchroMed II pump after an MRI to ensure continuation of therapy.
- Educate patients, caregivers, and family members to recognize the signs and symptoms associated with intrathecal drug therapy underdose or withdrawal. Patients receiving intrathecal baclofen therapy (e.g., Lioresal Intrathecal) are at higher risk for adverse events, as baclofen withdrawal can lead to a life-threatening condition if not treated promptly and effectively.

Customer Required Actions:

- Share this notice with all those who need to be aware of this issue within your organization and maintain a copy of this notice in your records.
- Please complete and return the Customer Acknowledgment Form enclosed with this letter, acknowledging that you have received this information.

Additional Information:

Medtronic has notified the Competent Authority of your country of this action.

We regret any inconvenience this may cause. We are committed to patient safety and appreciate your prompt attention to this matter. If you have any questions regarding this communication, please contact your Medtronic Representative .

Sincerely,

Ayman Doughan
Business Manager

Enclosure(s):

- Appendix 1: Product Scope
- Customer Acknowledgment Form

Appendix 1: Product Scope

Product Number/ CFN	UDI-Device Identifier (GTIN/UPN)
PUMP 8637-20 SYNCHROMED II (8637-20)	00643169345201, 00643169384101, 00643169700925, 00643169700932, 00643169732278, 00643169999862, 00763000122676, 00763000421748, 00763000421755, 00763000689643
PUMP 8637-40 SYNCHROMED II (8637-40)	00643169345232, 00643169384170, 00643169384231, 00643169701021, 00643169701038, 00643169701090, 00643169732377, 00643169999961, 00763000122768, 00763000421830, 00763000421847, 00763000689575

FA1367 Customer Acknowledgement Form - Response is required SynchroMed™ II MRI Guideline Update

Please complete this Form in its entirety.

Date: _____

Name of Person Completing this Form: _____

Title: _____

Direct Phone #: _____

Email: _____

Account Name: _____

Account Number: _____

Account Address: _____

City: _____ Zip Code: _____

Country: _____

I have read and understand the instructions provided and acknowledge receipt of the notification regarding the use of the **SynchroMed™ II** by signing below. I also agree to further distribute and communicate this important information within my facility and to anyone whom I have further distributed **SynchroMed™ II** as required.

Name: (print) Signature: Date:

If you have any questions regarding this notification, please contact your Medtronic sales representative.

PLEASE EMAIL OR FAX THIS ACKNOWLEDGEMENT TO:

nahar.s.alsurayi@medtronic.com