



To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No 243 dated 29/12/2022 Regarding NCMDR recall of Mahurkar Acute Dual Lumen High Flow (13.5 French) Hemodialysis Catheters (Mahurkar QPlus) from (mfr: Covidien LLC).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information





Circular No. 243/2022


05-06-1444 H

29 -12-2022

نتقدم بثقة
Moving Forward
with Confidence



Recall of Dental products - dental material from Ivoclar Vivadent AG.

Source	NCMDR- National Center for Medical Devices Reporting- SFDA https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=6&rid=18383
Product	VivaPen Snap-On Cannula and Adhese Universal in the VivaPen.
Description	Dental products - dental material.
Manufacturer	Ivoclar Vivadent AG.
The affected products	Product Codes 752382, 752383, 752384, 745762, 745763 all batches
Reason	The affected cannulas may not be reliably fastened to the VivaPen after having attached them.
Action	<ol style="list-style-type: none">1. Attach the VivaPen cannula by snapping it into place as indicated in the Instructions for Use. Make sure that you tangibly overcome a point of resistance when attaching the cannula.2. If you do not feel a "snap" when attaching the cannula, do not use it and discard it. If there are concerns regarding the fit of the cannula, please contact your distributor rep and you will be provided with a free replacement as quickly as possible.3. Always use the VivaPen with the respective protective sleeve (see Instructions for Use, chapter 2.1.4).4. Correct use of the protective sleeve will prevent the cannula from coming off inside the patient's mouth and being inhaled or ingested.5. Sensitize the users/customers to consider the described risk-minimizing measures for all batches of VivaPen cannulas. Please contact the supplier for the necessary information.6. Return all unused affected items.7. Contact the local agent for remedial action.
Product image	<p>1. Attach the VivaPen cannula by snapping it into place as indicated in the Instructions for Use. Make sure that you tangibly overcome a point of resistance when attaching the cannula.</p> 
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: Med-device@moh.gov.om

Dr. Mohammed Hamdan Al Rubaie

Director General

