



رؤية عمان 2040  
نتقدم بثقة  
Moving Forward  
with Confidence

To:

**THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES**

**Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)**

**Director General of Engineering Affairs, MOH**

**Director General of Royal Hospital**

**Director General of Khoula Hospital**

**Director General of Medical Supplies (MOH)**

**Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)**

**Hospital Director (Al Nahda Hospital)**

**Hospital Director (Al Massara Hospital)**

**The Head of Medical Services in SQU Hospital**

**The Head of Medical Services in Royal Oman Police**

**The Head of Medical Services in Ministry of Defence**

**The Head of Medical Services in The Diwan**

**The Head of Medical Services in The Sultan's Special Force**

**The Head of Medical Services in Internal Security Services**

**The Head of Medical Services in Petroleum Development of Oman**

**The Head of Medical Services in LNG Oman**

**ALL PRIVATE PHARMACIES & DRUG STORES**

After Compliments,

Please find attached our Circular No 220 dated 29/11/2022 Regarding NCMDR Field Safety Corrective Action of V. Mueller Bipolar Forceps from (mfr: Becton Dickinson & Co. (BD)).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information



**PADC**  
المديرية العامة للصيدلة والرقابة الدوائية  
Directorate General of Pharmaceutical  
Affairs & Drug Control



ص.ب: 393 مسقط - الرمز البريدي: 100 - هاتف: 22357111 - فاكس: 22358489

P.O. Box: 393 Muscat - Postal Code: 100 - Tel: 22357111 - Fax: 22358489

dgpa\_dc Email: dg-padc@moh.gov.om



Circular No. 220/2022

5 -05-1444 H

29 -11-2022

نتقدم بثقة  
Moving Forward  
With Confidence

رؤية عمان 2040  
Vision Oman 2040

**Field Safety Corrective Action of V. Mueller Bipolar Forceps from Becton Dickinson & Co. (BD).**

|                       |  |
|-----------------------|--|
| Source                | NCMDR- National Center for Medical Devices Reporting- SFDA<br><a href="https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=17325">https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=17325</a>  |
| Product               | V. Mueller Bipolar Forceps.  |
| Description           | Bipolar Forceps.   |
| Manufacturer          | Becton Dickinson & Co. (BD).   |
| The affected products | REF: See "Table 1" in the attached FSN.<br>Lot Numbers: All lot numbers.   |
| Reason                | The Instructions For Use (IFU) provided by BD is missing certain content that is contained in the manufacturer's IFU. This missing content is related to the interface with the device power supply, as well as cleaning and maintenance instructions.                         |
| Action                | 1. Replace your current IFU for the products listed in Table 1 with the updated IFU provided with this notification.<br>2. Contact the local agent for remedial action.  |
| Comments              | Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: <a href="mailto:Med-device@moh.gov.om">Med-device@moh.gov.om</a> |

Dr. Mohammed Hamdan Al Rubaie  
Director General



**PADC**  
المديرية العامة للصيدلة والرقابة الدوائية  
Directorate General of Pharmaceutical  
Affairs & Drug Control



ص.ب: ٣٩٣ مسقط - الرمز البريدي: ١٠٠ - هاتف: ٢٢٣٥٧١١١ - فاكس: ٢٢٣٥٨٤٨٩

P.O. Box: 393 Muscat - Postal Code: 100 - Tel: 22357111 - Fax: 22358489

dgpa\_dc Email: dg-padc@moh.gov.om





BD Switzerland Sàrl  
Terre Bonne Park – A4  
Route de Crassier 17  
1262 Eysins – Switzerland  
Tél: +41 21 556 30 Fax:  
+41 21 556 30 99  
www.BD.com

7<sup>th</sup> November 2022

**PRODUCT NOTIFICATION – SUR-22-4498**

**V. Mueller™ Bipolar Forceps**

**REF:** See Table 1 **Lot Numbers:** All lot numbers

**Type of Action:** Advisory

Dear Customer,

BD is issuing an advisory product notification for all lot numbers of **V. Mueller™ Bipolar Forceps**. According to our distribution records your organisation may have received the impacted product in Table 1. Product was distributed between December 2011 and January 2022.

| Product Code (REF) | Product Name   | Lot Number | UDI            |
|--------------------|--|------------|----------------|
| F-1305             | V. Mueller™ SCOVILLE-GREENWOOD BAYONET IRRIGATING BIPOLAR FORCEPS INSULATED, 1.5MM TIP OVERALL LENGTH 7-3/4" (195MM) | All Lots   | 10885403040856 |
| F-5125             | V. Mueller™ TITANIUM BAYONET BIPOLAR FORCEPS, STRAIGHT, FINE 0.5MM TIP OVERALL LENGTH 8-3/4" (225MM)                 | All Lots   | 10885403041211 |
| F-5126             | V. Mueller™ TITANIUM BAYONET BIPOLAR FORCEPS, 1.0MM TIP, STRAIGHT MEDIUM OVERALL LENGTH 8-3/4" (225MM)               | All Lots   | 10885403041228 |
| F-1000             | V. Mueller™ CUSHING BAYONET BIPOLAR INSULATED FORCEPS, 0.7MM TIP OVERALL LENGTH 7-1/2" (190MM)                       | All Lots   | 10885403040702 |
| F-3015             | V. Mueller™ ADSON BIPOLAR FORCEPS 1.0MM TIP, INSULATED OVERALL LENGTH 4-3/4" (120MM)                                 | All Lots   | 10885403040979 |
| F-1035             | V. Mueller™ HARDY-STYLE BAYONET BIPOLAR FORCEPS, 0.5MM TIP, INSULATED OVERALL LENGTH 8-1/2" (215MM)                  | All Lots   | 10885403040740 |
| F-5022             | V. Mueller™ TITANIUM ROUND HANDLE BAYONET BIPOLAR FORCEPS 1.5MM TIP OVERALL LENGTH 10-1/2" (265MM)                   | All Lots   | 10885403041075 |
| F-5018             | V. Mueller™ TITANIUM ROUND HANDLE BAYONET BIPOLAR FORCEPS, 1.0MM TIP OVERALL LENGTH 10-1/2" (265MM)                  | All Lots   | 10885403041068 |
| F-5016             | V. Mueller™ TITANIUM ROUND HANDLE BAYONET BIPOLAR FORCEPS 0.7MM TIP OVERALL LENGTH 10-1/2" (265MM)                   | All Lots   | 10885403041051 |
| F-5067             | V. Mueller™ ROUND HANDLE BAYONET BIPOLAR FORCEPS 0.7MM TIP OVERALL LENGTH 10-1/2" (265MM)                            | All Lots   | 10885403041167 |



BD Switzerland Sàrl  
Terre Bonne Park – A4  
Route de Crassier 17  
1262 Eysins – Switzerland  
Tél: +41 21 556 30 Fax:  
+41 21 556 30 99  
www.BD.com

|        |   |          |                |
|--------|---|----------|----------------|
| F-5002 | V. Mueller™ TITANIUM ROUND HANDLE BAYONET BIPOLAR FORCEPS, 1.0MM TIP OVERALL LENGTH 8-1/2" (215MM)    | All Lots | 10885403040986 |
| F-5073 | V. Mueller™ ROUND HANDLE BAYONET BIPOLAR FORCEPS 1.5MM TIP OVERALL LENGTH 10-1/2" (265MM)             | All Lots | 10885403041198 |
| F-5127 | V. Mueller™ TITANIUM BAYONET BIPOLAR FORCEPS, 0.5MM TIP, ANGLED-UP FINE OVERALL LENGTH 8-3/4" (225MM) | All Lots | 10885403041235 |
| F-2002 | V. Mueller™ CUSHING BIPOLAR FORCEPS 1.5MM TIP, INSULATED OVERALL LENGTH 7" (180MM)                    | All Lots | 10885403040900 |

**Table 1: Impacted product**

This advisory is limited to the product codes listed in Table 1. No other product codes are affected.

### **Description of the Problem**

BD has identified that the Instructions For Use (IFU) provided by BD is missing certain content that is contained in the manufacturer's IFU. This missing content is related to the interface with the device power supply, as well as cleaning and maintenance instructions. Refer to the attached updated **V. Mueller™ Bipolar Forceps IFU**.

The potential consequences of inappropriate power supply to the **V. Mueller™ Bipolar Forceps** device include the inability to effectively induce hemostasis, on one extreme, or electrical burn injuries, on the other. Either of these complications could adversely affect the efficient and safe delivery of care to patients.

The potential consequences of inappropriate cleaning, maintenance, and storage of the **V. Mueller™ Bipolar Forceps** device include contamination of sterile operative fields leading to infection, abscess formation, sepsis and, at an extreme, septic shock.

To date there has been no adverse events related to this issue.

**There is no requirement for customers to return any V. Mueller™ Bipolar Forceps to BD. These products can continue to be used in accordance with the guidance in this notification.**

### **Actions taken by BD:**

BD will implement appropriate corrective actions to prevent recurrence and is providing, with this notification, the **V. Mueller™ Bipolar Forceps IFU** which contains the missing content.



BD Switzerland Sàrl  
Terre Bonne Park – A4  
Route de Crassier 17  
1262 Eysins – Switzerland  
Tél: +41 21 556 30 Fax:  
+41 21 556 30 99  
www.BD.com

#### **Customer Actions:**

- Replace your current IFU for the products listed in Table 1 with the updated IFU provided with this notification.
- Complete and return the Customer Response Form **even if you no longer have any inventory remaining in your facility by 30<sup>th</sup> November 2022.**
- Circulate this notification to all those who need to be aware within your organization or to any organization where the potentially affected devices have been transferred.
- If you experience any issues with the **V. Mueller™ Bipolar Forceps**, please report as a complaint as per your normal process.

#### **Distributor Actions:**

- Review the information in **Table 1** and determine if **V. Mueller™ Bipolar Forceps** in your possession are impacted.
- Replace the current IFU with the updated IFU provided with this product notification.
- Identify the facilities where you have distributed affected product and notify them immediately of this notification. Have your customers complete and return the Customer Response form to your organisation for reconciliation purposes by **30<sup>th</sup> November 2022.**
- Complete and return the Customer Response Form following completion of your reconciliation activities.

|  | <b>End User with Inventory</b>  | <b>End User with ZERO inventory</b>   | <b>Where to send completed form</b>                                |
|--|---|---|--|
| Purchased <b>directly</b> from BD                        | Complete the form in its entirety and ensure that all recommended actions have been implemented as required | Complete the form in its entirety and retain a copy of this notification for your records | <a href="mailto:EMEAFieldAction@bd.com">EMEAFieldAction@bd.com</a> |
| Purchased from a <b>distributor/3<sup>rd</sup> party</b> | Complete the form in its entirety   | Complete the form in its entirety and retain a copy of this notification for your records | Return the form to your distributor/3 <sup>rd</sup> party          |



BD Switzerland Sàrl  
Terre Bonne Park – A4  
Route de Crassier 17  
1262 Eysins – Switzerland  
Tél: +41 21 556 30 Fax:  
+41 21 556 30 99  
[www.BD.com](http://www.BD.com)

**Contact Reference Person**

BD is committed to advancing the world of health. Our primary objectives are patient safety and user safety and providing you with quality products. We apologise for the inconvenience this situation may cause you and thank you in advance for helping BD to resolve this matter as quickly and effectively as possible.

Sincerely,

A handwritten signature in cursive script, appearing to read "L. Darrock".

Lorna Darrock  
Associate Director, Post Market Quality  
EMA Quality