



To:

**THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES**

**Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)**

**Director General of Engineering Affairs, MOH**

**Director General of Royal Hospital**

**Director General of Khoula Hospital**

**Director General of Medical Supplies (MOH)**

**Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)**

**Hospital Director (Al Nahda Hospital)**

**Hospital Director (Al Massara Hospital)**

**The Head of Medical Services in SQU Hospital**

**The Head of Medical Services in Royal Oman Police**

**The Head of Medical Services in Ministry of Defence**

**The Head of Medical Services in The Diwan**

**The Head of Medical Services in The Sultan's Special Force**

**The Head of Medical Services in Internal Security Services**

**The Head of Medical Services in Petroleum Development of Oman**

**The Head of Medical Services in LNG Oman**

**ALL PRIVATE PHARMACIES & DRUG STORES**

After Compliments,

Please find attached our Circular No 171 dated 21/5/23 Regarding NCMDR  
FSCA of Portrait Wearable Pulse Oximetry Sensors from (mfr: GE Healthcare).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information





ننمذ بثقة  
Moving Forward  
with Confidence




Circular No. 171 / 2023

05-02-1445 H

21 -08-2023

Field Safety Corrective Action of Portrait Wearable Pulse Oximetry Sensors from GE Healthcare

Source	NCMDR- National Center for Medical Devices Reporting- SFDA <a href="https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=19650">https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=19650</a>
Product	Portrait Wearable Pulse Oximetry Sensors.
Description	Pulse Oximetry Sensors.
Manufacturer	GE Healthcare.
Local Agent	Muscat Pharmacy & Stores LLC.
The affected products	ITEM ( PRODUCT CODE; REF #; GTIN): Portrait SpO2 P-W01 (SRY; 2096441-001; 00195278264800) Portrait SpO2 P-SE01 (SRY; 2104044-001; 00195278264824)
Reason	If above affected devices are saturated with fluids, can lead to possible loss of SpO2 monitoring.
Action	<ol style="list-style-type: none"><li>1. Visually inspect the affected devices for damage. Damaged devices should not be used (see Figure 1 in the attachment).</li><li>2. If the 'Faulty SpO2 probe' alarm occurs, replace the sensor as described in Sections 7 and 18 of the Portrait™ Mobile Monitoring Solution User Manual and discard the failed sensor.</li><li>3. Avoid saturating the affected devices with fluids. For example, remove the sensor before contact with water (e.g. hand washing and taking a shower). If the sensor comes in contact with fluids, dry the sensor completely before continuing use.</li><li>4. Follow the Portrait™ Mobile Monitoring Solution User Manual instructions (Section 15) for cleaning, disinfection and care.</li><li>5. Contact the local agent for remedial action.</li></ol>
comments	 <p>Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: <a href="mailto:Med-device@moh.gov.om">Med-device@moh.gov.om</a></p>

Dr. Mohammed Hamdan Al Rubaie

Director General

ص.ب: 393 مسقط - الرمز البريدي: 100 - هاتف: 22357111 - فاكس: 22358489

P.O. Box: 393 Muscat - Postal Code: 100 - Tel: 22357111 - Fax: 22358489

dgpa\_dc Email: [dg-padc@moh.gov.om](mailto:dg-padc@moh.gov.om)