



نقدم بثقة  
Moving Forward  
with Confidence



To:

**THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES**  
**Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)**  
**Director General of Engineering Affairs, MOH**  
**Director General of Royal Hospital**  
**Director General of Khoula Hospital**  
**Director General of Medical Supplies (MOH)**  
**Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)**  
**Hospital Director (Al Nahda Hospital)**  
**Hospital Director (Al Massara Hospital)**  
**The Head of Medical Services in SQU Hospital**  
**The Head of Medical Services in Royal Oman Police**  
**The Head of Medical Services in Ministry of Defence**  
**The Head of Medical Services in The Diwan**  
**The Head of Medical Services in The Sultan's Special Force**  
**The Head of Medical Services in Internal Security Services**  
**The Head of Medical Services in Petroleum Development of Oman**  
**The Head of Medical Services in LNG Oman**  
**ALL PRIVATE PHARMACIES & DRUG STORES**

After Compliments,

Please find attached our Circular No 109 dated 28/7/2024 Regarding NCMDR Recall of BLUperc® Percutaneous Dilation Tracheostomy Kit from (mfr: Smith's Medical International Limited).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DSC
- Director of Pharmacovigilance & Drug Information Dept, DSC
- Director of Drug Control Department, DSC
- Director of Pharmaceutical Licensing Department, DSC
- Director of Central Quality Control Lab., DSC
- Supdt. of Central Drug Information



Circular No. 109/2024

22-01-1446 H  
28-07-2024

نتقدم بثقة  
Moving Forward  
with Confidence



Recall of BLUperc® Percutaneous Dilation Tracheostomy Kit from Smith's Medical International Limited.

Source	NCMDR - National Center Medical Device Reporting- SFDA. <a href="https://ncmdr.sfda.gov.sa/Secure/CA/CAViewRecall.aspx?caid=4&amp;rid=21122">https://ncmdr.sfda.gov.sa/Secure/CA/CAViewRecall.aspx?caid=4&amp;rid=21122</a>
Product	BLUperc® Percutaneous Dilation Tracheostomy Kit.
Description	Tracheotomy.
Manufacturer	Smith's Medical International Limited.
Local agent	Muscat Pharmacy & Stores LLC.
The affected products	SKU/List/Model Number (Lot/Serial Number): 101/561/090 (4439066)
Reason	The potential for a disconnection of the pilot balloon from the tracheostomy inflation line within specific lots.
Action	1. Please discontinue use and discard all stock of the affected lot. If, due to extenuating circumstances, an alternative device cannot be used, all instructions, including warnings and cautions contained in the Instructions for Use Documentation must be followed with heightened awareness when checking the device prior to use and while in use. 2. Contact the local agent for remedial action.
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: <a href="mailto:vigilance-md@moh.gov.om">vigilance-md@moh.gov.om</a>

Dr. Mohammed Hamdan Al Rubaie  
Director General



# URGENT: FIELD SAFETY NOTICE

**BLUselect® Tracheostomy Tube Kits , BLUselect® Suctionaid® Tracheostomy Tube Kits, BLUgriggs® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with BLUselect® Tracheostomy Tube with or without Forceps, BLUperc® Dilation Procedural Tray with Single Stage Dilator Products, BLUperc® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with or without BLUselect® Tracheostomy Tube**


04 July 2024

Dear Valued BLUSelect® Customers,

Smiths Medical, Inc. is issuing this Urgent Field Safety Notice to notify you of a potential defect with the following BLUSelect®, BLUgriggs® and BLUperc® products listed in *Attachment 1\_Affected Product*. This letter details the issue and the required steps for you to complete.

**Issue:**

Smiths Medical has identified the potential for a disconnection of the pilot balloon from the tracheostomy inflation line within specific lots of the BLUSelect®, BLUgriggs® and BLUperc® products because of a manufacturing defect. See the photo example of the issue below.



**Potential Risk:**

If the pilot balloon used to inflate the tracheostomy cuff becomes detached from the inflation line, the cuff pressure may not be maintained which can lead to inadequate ventilation of, and increased risk of aspiration to the patient. To date, Smiths Medical has received thirteen (13) reports of serious injury and zero (0) deaths associated with this issue.

**Affected Product:**

SKU/List/Model Number	Lot/Serial Number	Shipped Quantity
101/561/090	4439066	1

**Smiths Medical Actions:**

Smiths Medical is sending this notification to all BLUSelect®, BLUgriggs® and BLUperc® customers who received products from Smiths Medical listed in *Attachment 1\_Affected Product*. Smiths Medical will provide replacement product(s) or credit to affected customers upon receipt of a completed response form to certify product destruction.

**Customer Required Actions:**

When using the device, all instructions, including warning and cautions contained in the Instructions for Use Documentation must be followed with heightened awareness. Please complete the following actions listed below

1. Check all inventory locations within your institution for the affected catalog numbers and lot numbers listed in the notification and discontinue use. Discard all affected products following your institution's process for discarding. If discarding is not immediately possible at your facility, then the product should be quarantined until disposal.
2. Share this notification with all potential users of the device to ensure they are aware of this notification and proposed mitigations. If the devices are used at another location, please ensure this communication is delivered there.
3. Complete and return the attached Customer Response Form to [EMEA-FSN@icumed.com](mailto:EMEA-FSN@icumed.com) within 10 days of receipt to acknowledge your understanding of this notification.
4. **DISTRIBUTORS:** If you have distributed potentially affected products to your customers, please immediately forward this notice to them and request that they complete the response form and return it to **YOU**. Then the **DISTRIBUTOR** must complete a SINGLE form with the required details and return to [EMEA-FSN@icumed.com](mailto:EMEA-FSN@icumed.com)

For further inquiries, please contact Smiths Medical using the following information:

Smiths Medical Contact	Contact Information	Areas of Support
Global Complaint Management	<a href="mailto:globalcomplaints@icumed.com">globalcomplaints@icumed.com</a>	To report adverse events or product complaints
Customer Service	<a href="https://www.icumed.com/about-us/contact-us">https://www.icumed.com/about-us/contact-us</a>	Questions about product replacement and/or credit.

**Your country regulatory agency has been notified of this action**

Smiths Medical is committed to patient safety and is focused on providing exceptional product reliability and the highest level of customer satisfaction. Thank you for your prompt support on this important matter. We appreciate your cooperation.

Sincerely,



Andy Mathein  
Vice President of Quality

**Enclosures:**

- Customer Response Form (See Below)
- Affected Product List (Attachment 1)

# URGENT: FIELD SAFETY NOTICE – RESPONSE FORM

**BLUselect® Tracheostomy Tube Kits , BLUselect® Suctionaid® Tracheostomy Tube Kits, BLUgriggs® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with BLUselect® Tracheostomy Tube with or without Forceps, BLUperc® Dilation Procedural Tray with Single Stage Dilator Products, BLUperc® Percutaneous Dilation Tracheostomy Procedural Kit or Tray with or without BLUselect® Tracheostomy Tube**

4 July 2024

Check your inventory and complete the information below, even if you do not have the affected product. Failure to complete all sections of this page may result in improper, delayed or denied credit.

Please return the completed form to [EMEA-FSN@icumed.com](mailto:EMEA-FSN@icumed.com), If you have questions about this form please contact [EMEA-FSN@icumed.com](mailto:EMEA-FSN@icumed.com) or your local sales representative.

Name of Hospital / Facility	
Hospital / Facility Address	
Telephone Number	
Name and Title of Person Completing this Form	
Signature of Person Completing this Form	
Date	
If Purchased through a distributor, please list distributor name/location here for traceability purposes	

Please select one:

- I have **NO** affected products (complete and return this form to the e-mail address above)
- YES**, I have affected products, I have notified users in my facility and I have followed the instructions provided to me and destroyed all affected items (see table below)

If you have affected product on hand, please complete table 1 below:

**TABLE 1**

Item / SKU Number	Lot Number	Quantity in inventory (Eaches)	Quantity Destroyed (Eaches)	Date of Destruction

If you have distributed the product further, please complete table 2 below with collated information received from your customers and respond to ICU Medical with the overall information.

**TABLE 2**

Item / SKU Number	Lot Number	Quantity destroyed locally (Eaches)	Date of Destruction

Adverse events and complaints associated with the use of this product should be reported and emailed to Smiths Medical's Global Complaint Management Department at [globalcomplaints@icumed.com](mailto:globalcomplaints@icumed.com).

## ADDITIONAL AFFECTED PRODUCT DESTROYED

Item / SKU Number	Lot Number	Quantity in inventory (Eaches)	Quantity Destroyed (Eaches)	Date of Destruction

## Medical Devices Sector

قطاع الأجهزة الطبية

- Home
- Published FSNs/Recalls
- About NCMDR
- Contact Us
- FAQ
- Login

# NCMDR

National Center for Medical Devices Reporting


المركز الوطني لبلاغات الأجهزة والمنتجات الطبية

## NCMDR Recall

**Reference Number:** mdprc 016 07 24 000

[Back](#)

**Date submitted:** 7/9/2024

<b>Manufacturer:</b>	Smith's Medical International Limited
<b>Device Type:</b>	BLUperc® Percutaneous Dilation Tracheostomy Kit
<b>Description:</b>	Tracheotomy
<b>Medical Device Identifier:</b>	SKU/List/Model Number (Lot/Serial Number): 101/561/090 (4439066)
<b>Reason of Field Safety Corrective Action:</b>	The potential for a disconnection of the pilot balloon from the tracheostomy inflation line within specific lots.
<b>Remedy Action:</b>	Please discontinue use and discard all stock of the affected lot. If, due to extenuating circumstances, an alternative device cannot be used, all instructions, including warnings and cautions contained in the Instructions for Use Documentation must be followed with heightened awareness when checking the device prior to use and while in use.
<b>Athorized Representative/Importer/Distributor:</b>	Almadar medical Est.
<b>Report Source:</b>	NCMDR
<b>Source Ref. Number:</b>	SA-08-07-24-499
<b>SFDA Comments:</b>	SFDA urges all healthcare providers that have devices subjected to this safety alert to contact the company.
<b>Attachments:</b>	 <a href="#">Smiths Medical International Limited.pdf</a>

[View History](#)

Copyright © 2008 Saudi Food and Drug Authority. All rights reserved.