

To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

ALL PRIVATE PHARMACIES & DRUG STORES


After Compliments,

Please find attached our Circular No 250 dated 29/12/2022 Regarding NCMDR Recall of NEO-fit Neonatal Endotracheal Tube Grip of MR systems from (mfr: Cooper Surgical Inc).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information





Circular No. 2501/2022


05 -06-1444 H

29 -12-2022

نتقدم بثقة
Moving Forward
with Confidence



Recall of NEO-fit Neonatal Endotracheal Tube Grip from Cooper Surgical Inc.

Source	NCMDR- National Center for Medical Devices Reporting- SFDA https://ncmdr.sfda.gov.sa/Secure/CA/CAViewRecall.aspx?caid=6&rid=18394
Product	NEO-fit Neonatal Endotracheal Tube Grip.
Description	Neonatal Endotracheal Tube Grip.
Manufacturer	Cooper Surgical Inc.
Local agent	Muscat Pharmacy & Stores LLC.
The affected products	Part Number: 42-2540 Lots: Refer to "Appendix A" in the attached FSN.
Reason	Complaints relating to the metal clip becoming detached from the product strap and there is the risk that this clip may be ingested by the infant due to its proximity to the oral cavity.
Action	1. Return affected unexpired product. 2. Contact the local agent for remedial action.
Product image	
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: Med-device@moh.gov.om

Dr. Mohammed Hamdan Al-Rubaie

Director General





75 Corporate Drive
Trumbull, CT 06611

T 203 601 5200
www.coopersurgical.com

CooperSurgical®

Bundesinstitut für Arzneimittel und Medizinprodukte (BfArM)
Kurt-Georg-Kiesinger-Allee 3
53175 Bonn
Germany
Phone: +49 (0)228-207-30
Email: poststelle@bfarm.de

10/18/2022

Medical Device Field Safety Corrective Action:
NEO-fit™ Neonatal Endotracheal Tube Grip
Part Number: 42-2540
Lots: Appendix A

Dear Sir/Madam

CooperSurgical will send a Field Safety Notice (FSN) notification to all affected customers. The products are being recalled due to receiving complaints relating to the metal clip becoming detached from the product strap and there is the risk that this clip may be ingested by the infant due to its proximity to the oral cavity.

CooperSurgical will notify affected consignees via certified mail with confirmed delivery receipt. The letter will identify the affected product and detail the necessary actions to be taken by customers who have received affected products. The letter will ask customers to contact CooperSurgical if they have any questions regarding the recall. At CooperSurgical's expense arrangements will be made to return affected unexpired product to CSI for credit.

There are up to 220,500 unexpired devices that are considered affected and distributed globally. CooperSurgical determined that the affected products were distributed to 2 consignees in Germany.

CooperSurgical has initiate CAPA 793 to determine root cause and implement a Corrective Action Plan (CAP) to address the root cause and prevent recurrence. Upon completion of all recall activities, CooperSurgical will formally submit a closure report. Copies of the Health Hazard Evaluation for risk analysis is enclosed.

Sincerely,

Manager of Risk Management

Appendix A – Lot Numbers of Potentially Affected Devices (Part Number: 42-2540)

Lot Number	Lot Number	Lot Number	Lot Number
247020	292892	298064	307729
269476	292893	298065	307730
269477	292894	301600	308590
269478	292895	301601	308591
269479	292896	301602	309456
269480	292897	301603	309544
269481	292898	301604	309910
269482	292899	301605	309912
269483	292900	301606	309913
272992	292901	301607	309914
274957	292902	302944	311350
274958	292903	303769	311465
274959	292904	303771	311533
274960	296234	303772	311690
274961	296235	303949	312729
274962	296236	303950	312730
274963	296237	304053	312731
286859	296238	304363	313142
286861	296239	304364	313273
286862	296240	305016	313274
286863	296241	305017	314493
291949	296890	305018	314675
291950	296891	305019	314676
291951	296892	305020	314969
291952	296893	305021	315504
291953	297033	306282	315505
292878	297034	306339	315626
292879	297035	306340	315810
292880	297036	306341	316162
292881	297037	306342	316866
292882	298047	306343	316976
292883	298048	307080	317032
292884	298049	307081	318970
292885	298051	307616	319205
292886	298053	307723	319298
292887	298054	307724	320400
292888	298055	307725	321242
292889	298056	307726	
292890	298062	307727	
292891	298063	307728	