



To:

**THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES**  
Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)  
Director General of Engineering Affairs, MOH  
Director General of Royal Hospital  
Director General of Khoula Hospital  
Director General of Medical Supplies (MOH)  
Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)  
Hospital Director (Al Nahda Hospital)  
Hospital Director (Al Massara Hospital)  
The Head of Medical Services in SQU Hospital  
The Head of Medical Services in Royal Oman Police  
The Head of Medical Services in Ministry of Defence  
The Head of Medical Services in The Diwan  
The Head of Medical Services in The Sultan's Special Force  
The Head of Medical Services in Internal Security Services  
The Head of Medical Services in Petroleum Development of Oman  
The Head of Medical Services in LNG Oman  
**ALL PRIVATE PHARMACIES & DRUG STORES**

After Compliments,

Please find attached our Circular No 247 dated 29/12/2022 Regarding NCMDR Recall of ID-DiaPanel from (mfr: DiaMed GmbH).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information





Circular No. 247/2022

نحن نقدم  
مضيافاً  
تقدمنا  
بثقة  
Moving Forward  
With Confidence



05-06-1444 H

29-12-2022

### Recall of ID-DiaPanel from DiaMed GmbH.

|                       |  |
|-----------------------|--|
| Source                | NCMDR - National Center Medical Device Reporting- SFDA.<br><a href="https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=10&amp;.rid=18374">https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=10&amp;.rid=18374</a>   |
| Product               | ID-DiaPanel.   |
| Description           | IVD.   |
| Manufacturer          | DiaMed GmbH.   |
| Local agent           | Bahwan Healthcare Center.  |
| The affected products | ID-DiaPanel (1-11)<br>Catalogue number: 004114<br>Lot numbers: 769898541 & 45161541<br>Expiry date: 26/12/2022   |
| Reason                | Inability of the IH-500 instrument to read the 2D barcode of the ID-DiaPanel when performing Antibody Identification tests for patients and/or donors on the IH-500 instrument.  |
| Action                | 1. Refer to "Recall Action Instructions" in the attached recall.<br>2. Contact the local agent for remedial action.  |
| comments              | Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: <a href="mailto:Med-device@moh.gov.om">Med-device@moh.gov.om</a> |

Dr. Mohammed Hamdan Al Rubaie

Director General



**PADDC**  
المديرية العامة للصيدلة والرقابة الدوائية  
Directorate General of Pharmaceutical  
Affairs & Drug Control



ص.ب: 393 مسقط - الرمز البريدي: 100 - هاتف: 22357111 - فاكس: 22358489

P.O. Box: 393 Muscat - Postal Code: 100 - Tel: 22357111 - Fax: 22358489

dgpa\_dc Email: [dg-padc@moh.gov.om](mailto:dg-padc@moh.gov.om)



**Australian Government**

**Department of Health**

Therapeutic Goods Administration

# Recall Action Notification

ID-DiaPanel (1-11). An in vitro diagnostic medical device  
(IVD)

© Commonwealth of Australia 2022.

This work is copyright. You may reproduce the whole or part of this work in unaltered form for your own personal use or, if you are part of an organisation, for internal use within your organisation, but only if you or your organisation do not use the reproduction for any commercial purpose and retain this copyright notice and all disclaimer notices as part of that reproduction. Apart from rights to use as permitted by the Copyright Act 1968 or allowed by this copyright notice, all other rights are reserved and you are not allowed to reproduce the whole or any part of this work in any way (electronic or otherwise) without first being given specific written permission from the Commonwealth to do so. Requests and inquiries concerning reproduction and rights are to be sent to the TGA Copyright Officer, Therapeutic

## Important information on the System for Australian Recall Actions

The TGA publishes information about therapeutic goods supplied in the Australian market that have been subject to a recall action in a publicly searchable database.

Recall action means action taken by the responsible entity (being the person who is responsible for taking the recall action) to resolve a problem with therapeutic goods supplied in the Australian market that have, or may potentially have, deficiencies relating to safety, quality, efficacy (performance) or presentation.

- Recall actions include: the permanent removal of therapeutic goods from supply in the market, the taking of corrective action in relation to therapeutic goods (such as repair, modification, adjustment or relabelling) and, in the case of medical devices that have been implanted into patients, the issuing of a hazard alert containing information for health practitioners on how to manage patients.
- More information about Australian recall actions is available at <http://tga.gov.au/safety/recalls-about.htm>
- If you are taking a medicine, using a medical device or have had a medical device implanted into you, that is the subject to a recall action, and you have any concerns you should seek advice from a health professional. <http://www.healthdirect.org.au/>

### About the release of this information

While reasonable care is taken to ensure that the information is an accurate record of recall actions that responsible entities have reported to the TGA or of which the TGA has become aware, the TGA does not guarantee or warrant the accuracy, reliability, completeness or currency of the information or its usefulness in achieving any purpose.

To the fullest extent permitted by law, including but not limited to section 61A of the Therapeutic Goods Act 1989, the TGA will not be liable for any loss, damage, cost or expense incurred in or arising by reason of any person relying on this information.

The information contained in the SARA database is released under s 61(5C) of the Therapeutic Goods Act 1989. Copyright restrictions apply to the System of Australian Recall actions (SARA) <http://tga.gov.au/about/website-copyright.htm>.

## Recall detail

|   |  |
|---|--|
| <b>Type of Product<sup>i</sup></b>                  | Medical Device   |
| <b>TGA Recall Reference<sup>ii</sup></b>            | RC-2022-RN-01445-1   |
| <b>Product Name/Description<sup>iii</sup></b>       | ID-DiaPanel (1-11). An in vitro diagnostic medical device (IVD)<br><br>Catalogue number: 004114<br><br>Lot numbers: 769898541 & 45161541<br><br>Expiry date: 26/12/2022<br><br>ARTG 213161<br>(Bio-Rad Laboratories Pty Ltd - Reagent red blood cell IVDs)   |
| <b>Recall Action Level<sup>iv</sup></b>             | Hospital   |
| <b>Recall Action Classification<sup>v</sup></b>     | Class II   |
| <b>Recall Action Commencement Date<sup>vi</sup></b> | 24/11/2022   |
| <b>Responsible Entity<sup>vii</sup></b>             | Bio-Rad Laboratories Pty Ltd   |
| <b>Reason / Issue<sup>viii</sup></b>                | <p>Bio-Rad Laboratories has identified an issue regarding inability of the IH-500 instrument to read the 2D barcode of the ID-DiaPanel when performing Antibody Identification tests for patients and/or donors on the IH-500 instrument. This defect results in the following error message for IH-500 users only: "Unreadable reagent barcode". This leads to the impossibility to load the reagent vial using the 2D barcode.</p> <p>Out of the set of 11 vials that constitute the reference 004114, only vial 11 is impacted. This problem does not affect the instrument itself.</p> |
| <b>Recall Action<sup>ix</sup></b>                   | Product Defect Correction  |

|   |   |
|---|---|
| <b>Recall Action Instructions<sup>x</sup></b> | <p>Bio-Rad are advising customer that under following instructions the product remains usable:</p> <ol style="list-style-type: none"> <li>1. The deactivation of the On Board Time function on IH-500. The instrument would use the valid 1D barcode available on the product label. For the details of the procedure (deactivation of On Board Time function on IH-500) please contact your local Bio-Rad techsupport contact. Please note that an advanced or Admin user level could be required. This procedure will lead to the loss of the On Board Time Management information (which defines how long the vial may be used on board), an option to manage the Shelf-Life of the loaded vial. The control of the expiry date of the reagent is not affected.</li> <li>2. Contact your local technical support to determine the appropriate solution (such as reproducing the 2D barcode).</li> </ol> <p>Bio-rad are also advising customers to use the next batch available from November 25th 2022 incase the device is not usable under any conditions describe above</p> |
| <b>Contact Information<sup>vi</sup></b>       | 1800 224 354 - Technical support - Biorad   |

## Footnotes

<sup>i</sup> Type of Product: Medicine, Medical Device, or Biological

<sup>ii</sup> TGA Recall Reference: Unique number given by the TGA

<sup>iii</sup> Product Name/Description: Brand name (including active ingredient for medicines) and may include generic reference for the kind of medical devices. Includes all necessary information such as affected: catalogue / model and / or batch / serial numbers.

<sup>iv</sup> Recall Action Level: The level to which the recall action is to be undertaken. This is based on the significance of the risk and the channels through which the goods have been distributed. The recall action levels are / Wholesale / Hospital / Retail / Consumer.

- **Wholesale** - includes wholesalers and state purchasing authorities.
- **Hospital** - includes nursing homes and institutions, hospital pharmacists, ambulance services, blood and tissue banks and laboratories as well as wholesale as appropriate.
- **Retail** - includes retail pharmacists, medical, dental and other health care professionals as well as wholesale and hospital as appropriate.
- **Consumer** - includes patients and consumers, as well as wholesale, hospital and retail levels as appropriate.

<sup>v</sup> Recall Action Classification<sup>\*\*</sup>: Recall actions of therapeutic goods are classified based on the potential risk the deficiency poses to patients / consumers. They are classified as Class I, Class II or Class III.

- **Class I** - A situation in which there is a reasonable probability that the use of, or exposure to, the deficient therapeutic good(s) will cause serious adverse health consequences or death.
- **Class II** - A situation in which use of, or exposure to, the deficient therapeutic good(s) may cause temporary or medically reversible adverse health consequences, or where the probability of serious adverse health consequences is remote.
- **Class III**- A situation in which use of, or exposure to, the deficient therapeutic good(s) is not likely to cause adverse health consequences.

<sup>vi</sup> Recall Action Commencement Date: The date the recall strategy and communication was agreed by the TGA.

<sup>vii</sup> Responsible Entity: Sponsor / Supplier / Importer responsible for the recall actions.

<sup>viii</sup> Reason / Issue: Reason for the recall action.

<sup>ix</sup> Recall Action<sup>\*\*</sup>: Recall action is an action taken to resolve a problem with a therapeutic good already supplied in the market for which there are issues or deficiencies in relation to safety, quality, efficacy (performance) or presentation.

There are four distinct recall actions – recall, product defect correction, hazard alert and product defect alert.

- **Recall** - The permanent removal of an affected therapeutic good from supply or use in the market.
- **Product defect correction** - Repair, modification, adjustment or re-labelling of a therapeutic good. The corrective action may take place at the user's premises or any other agreed location.
- **Hazard alert** - Information issued to healthcare professionals about issues or deficiencies relating to an implanted medical device or biological product and advice about the ongoing management of patients.
- **Product defect alert** - Information issued to raise awareness about issues or deficiencies for a therapeutic good where a recall action will result in interruption of patient treatment or a medicine shortage, including advice to reduce potential risks of using affected goods.

<sup>x</sup> Recall Action Instructions: What customers with affected goods should do.

<sup>xi</sup> Contact Information: Who the customer should contact for additional information and clarification regarding the recall action.

<sup>\*\*</sup> These definitions are applicable to the 2017 URPTG (Implemented from Jan 15 2018). Recall Action types and Recall Action Classifications prior to 15 Jan 2018 can be found at:

<https://www.tga.gov.au/sites/default/files/recalls-urptg-170412.pdf>