





## Ministry of Health

<b>Document Title</b>	Urgent and out of hours sample SOP
<b>Document Type</b>	Procedure
<b>Directorate/Institution</b>	Diagnostic Laboratories Services at Directorate General of Specialized Medical Care
<b>Targeted Group</b>	Medical laboratories
<b>Document Author</b>	Dr. Hanaa Al-raimi
<b>Designation</b>	Consultant medical microbiologist
<b>Document Reviewer</b>	Microbiology documents development team
<b>Designation</b>	Microbiology documents development team
<b>Release Date</b>	May 2023
<b>Review Frequency</b>	Three Years

<b>Validated by</b>		<b>Approved by</b>	
Name	Dr. Muna Habib	Name	Dr. Badryah Al Rashidi
Designation	Director Department Development & Control ( DGQAC)	Designation	Director General of Primary Health Care
Signature		Signature	
Date	May 2023	Date	June 2023

## **Contents Table:**

Acronyms:.....	4
.1 Purpose.....	5
2. Scope.....	5
3. Definitions.....	5
4. Procedure.....	5
5. Responsibilities.....	7
6. Document History and Version Control.....	9
7. References:.....	10

## Acknowledgment

The diagnostic laboratories services at the Directorate General of Specialized Medical Care (DGSMC) at Ministry of Health (MOH) would like to thank and appreciate the great effort of the Microbiology documents development team. Participated and contributed personnel are:

Member name	Institution	Designation
Dr.Mahmoud Al Subhi	Rustaq Hospital	Team Leader Consultant medical microbiologist
Ms. Zainab Al Hadhrami	Directorate General of Specialized Medical Care	Team Coordinator Senior technologist specialist A
Ms. Saleh Al Shukairi	Ibra Hospital	Senior technologist specialist A
Dr. Hanaa Al Auraimi	Royal Police of Oman Hospital	Consultant medical microbiologist
Dr. Nawal AL Kindi	Khoula Hospital	Consultant medical microbiologist
Dr. Al Warith Al Kharusi	Nizwa Hospital	Consultant medical microbiologist
Dr. Abdulrahman Al Mahrouqi	Ibri Hospital	Specialist microbiologist pathologist
Dr. Nada Al Tamimi	Al Massara Hospital	Consultant medical microbiologist
Dr. Wafaa Al Tamtami	Armed Forces Hospital	Consultant medical microbiologist

**Acronyms:**

BA	Blood agar
CA	Chocolate agar
MAC	MaConkey
ATCC	American Type Culture Collection
H&S	Health and Safety
ID	Identification
IQC	Internal Quality Control
MDRO	Multidrug Resistant Organism
MRSA	Methicillin Resistant Staph. Aureus
SOP	Standard operating procedure
TAT	Turnaround time
WHO	World Health Organization

## **1. Purpose**

The purpose of this procedure is to provide instructions on how to deal with urgent samples received in microbiology lab.

## **2. Scope**

This document is applicable for all medical laboratories under MOH and other collaborative governmental and non-governmental health institutions.

## **3. Definitions**

3.1 Urgent specimens: samples that need to be processed within 1 - 2 hours after sample collection.

## **4. Procedure**

### **4.1. Clinical background:**

Urgent specimens are specimen where results may change patient management (e.g meningitis or Necrotizing fasciitis). Urgent specimens need to be processed within 1 - 2 hours after collection. Result should be available within one to 2 hours and communicated to the ward as soon as is available.

Specimens of high importance to management where an early processing will allow an earlier result to be acted upon or precious unrepeatable specimens the quality of which will diminish too much if left for any longer.

These tests need to be available 24 hours a day and should be processed out of hours as priority.

### **4.2. Principle:**

Most specimens deteriorate if there is a significant delay between sample collection from the patient & processing. The delay will adversely affect test result (e.g. cell count).

### **4.3. Pre – analytical stage:**

#### **4.3.1. Sample**

4.3.1.1 The following specimens are processed as urgent ONCE received in the laboratory and result reported immediately:

- Cell count and Direct Smear for Sterile sites (CSF, pleural fluid, pericardial fluid, joint aspirate, etc.)

- Cerebrospinal fluids (CSFs): microscopy and set up for culture, csf cryptococcal antigen
- Cell count for urine
- Direct smear of positive blood culture
- Other requests for STAT processing of specimens must be approved by a Microbiologist

4.3.1.2 The following specimen must be processed out-of-hours as urgent testing

- Cerebrospinal fluids (CSFs): microscopy and set up for culture  
CSF: Cryptococcal antigen  
Urine microscopy and set up for culture (pediatric)
- Pus samples (not swabs) for microscopy
- Tissues / fluids / aspirates (including ascetic, joint, and CAPD fluid microscopy)
- Corneal scrape
- Throat swabs from cases of possible meningococcal disease
- All specimens from ICU & SCBU wards
- Blood cultures
- Urine pregnancy test from suspected ectopic pregnancy or patient for X-ray
- Catheters
- Deep organ abscesses

4.3.2. Material: NA

4.3.3. Safety precaution:

- All specimens need to be treated as potentially infectious. Standard procedures for handling of biohazard material must be followed at all times. Universal Precautions must be practiced at all stages of these procedures.

4.3.4. Quality control:

- Check the expiry dates of all media, reagents and stains before use.
- All media, reagents, kits, and stains **MUST** be quality controlled before use.
- Identification tests should be run with appropriate controls.

- Record the quality control results in the appropriate QC sheet.

#### 4.4. Analytical stage:

procedure for handling urgent specimens:

- Urgent sample can be ordered 24/7 hours including Holidays
- The Microbiology Lab should be notified by phone about the impending arrival of an urgent specimen.
- Specimens should be hand delivered to the microbiology lab reception.
- When the urgent specimen does not arrive at the Lab within reasonable time, the technologist should track down the specimen by calling the ward to find out where and when the specimen was sent.
- Once received to should be processed as urgent.

#### 4.1.2 Out of hours' samples

Samples listed above should be processed out of hours, other Specimens will be stored in the lab and processed the next morning.

#### 4.5. Post – analytical stage:

4.5.1. Reporting: Phone all result of urgent sample and release it in LIS system with proper documentation

#### 4.5.2. documentation of urgent specimen

All the following information should be documented on the Urgent Specimen Log sheet and released in LIS system:

Patient details	Date	Test /specimen type	Time requested	Time received	Result	Result informed by	Result informed to	Signature

## 5. Responsibilities

### 5.1. Laboratory management shall:

- To ensure the implementation of communication of critical results

- To coordinate with end-users for approval of the critical list and its related channels
- To collaborate with concerned department towards achievement of closed loop communication

5.2. Responsible staff shall:

- To ensure the adherence critical result communication procedure
- To facilitate the alternative channels once needed

5.3. Quality manager /officer

- To follow up the implementation of the procedure
- To monitor regularly communication of critical results and raise non-conformance with corrective action once needed.

5.4. All lab staff shall:

- To adhere to the procedure.
- To document communication as recommended
- To report communication failures or incident



## 6. Document History and Version Control

Version	Description	Review Date
1	Initial Release	May 2026

## 7. References

<b>Title of book/ journal/ articles/ Website</b>	<b>Author</b>	<b>Year of publication</b>	<b>Page</b>
Policy # QPCMI03001, STAT Specimens Procedure	University health network/mount Sinai hospital, department of microbiology		
Specimen Timeliness	Betsi cad walader University health board (		