





Ministry of Health

Document Title	Bile investigation SOP
Document Type	Procedure
Directorate/Institution	Directorate General of Specialized Medical Care
Targeted Group	Medical laboratories
Document Author	Dr Wafa Nasser Al Tamtami
Designation	Senior Specialist Medical Microbiology
Document Reviewer	Microbiology documents development team
Designation	Microbiology documents development team
Release Date	May 2023
Review Frequency	Three Year

Validated by		Approved by	
Name	Dr. Muna Habib	Name	Dr.Badryah Al Rashidi
Designation	Director Department Development & Control (DGQAC)	Designation	Director General of Primary Health Care
Signature		Signature	
Date	May 2023	Date	June 2023

Contents Table:

Acronyms:.....	4
1. Purpose.....	5
2. Scope.....	5
3. Definitions.....	5
4. Procedure.....	5
5. Responsibilities.....	11
6. Document History and Version Control.....	12
7. References:.....	13

Acknowledgment

The diagnostic laboratories services at the Directorate General of Specialized Medical Care (DGSMC) at Ministry of Health (MOH) would like to thank and appreciate the great effort of the Microbiology documents development team. Participated and contributed personnel are:

Member name	Institution	Designation
Dr.Mahmoud Al Subhi	Rustaque Hospital	Team Leader Consultant medical microbiologist
Ms. Zainab Al Hadhrami	Directorate General of Specialized Medical Care	Team Coordinator Senior technologist specialist A
Ms. Saleh Al Shukairi	Ibra Hospital	Senior technologist specialist A
Dr. Hanaa Al Auraimi	Royal Police of Oman Hospital	Consultant medical microbiologist
Dr. Nawal AL Kindi	Khoula Hospital	Consultant medical microbiologist
Dr. Al Warith Al Kharousi	Nizwa Hospital	Consultant medical microbiologist
Dr. Abdulrahman Al Mahrouqi	Ibri Hospital	Specialist microbiologist pathologist
Dr. Nada Al Tamimi	Al Massara Hospital	Consultant medical microbiologist
Dr. Wafaa Al Tamtami	Armed Forces Hospital	Senior specialist microbiologist

Contributors:

- Ms. Budoor AÁmri, Armed Forces Hospital

Acronyms:

AST	Antimicrobial Susceptibility Testing
ATCC	American Type Culture Collection
BA	Blood agar
BHS	Beta hemolytic Streptococci
BSC II	Biosafety Cabinet Class two.
BSL	Biosafety Level
BMS	Biomedical Scientist
CA	Chocolate agar
CO2	Carbon Dioxide
CoNS	Coagulase Negative Staphylococcus
GNB	Gram Negative Bacilli
H&S	Health and Safety
ID	Identification
IQC	Internal Quality Control
LIS	Laboratory Information System
MAC	MaConkey
MDRO	Multidrug Resistant Organism
MRSA	Methicilin Resistant Staph. Aureus
MTZ	Metronidazole disk
O2	Oxygen
SOP	Standard operating procedure
XLD	Xylose Lysine Deoxycholate Agar

1. Purpose

This SOP describes the methods of processing bile for bacteriological culture.

2. Scope

This document is applicable for all medical laboratories under MOH and other collaborative governmental and non-governmental health institutions.

3. Definitions

3.1 **Cholangitis** is the inflammation of the biliary ducts. It may present as ascending or suppurative cholangitis.

3.2 **Ascending cholangitis:** A partial obstruction of the biliary ducts combined with bacterial infection of the bile. Bacteria are intermittently shed into the bloodstream. Ascending cholangitis is a common cause of sepsis following liver transplantation.

3.3 **Suppurative cholangitis:** The complete obstruction of the infected biliary system. Increased biliary pressure causes constant shed of bacteria into the bloodstream.

3.4 **Cholecystitis** is the inflammation of the gall bladder that is usually due to an infection. The infection is often secondary to gallstones.

3.5 **Emphysematous cholecystitis** : is an acute infective cholecystitis that involves gas-forming organisms, most commonly *Clostridium perfringens*.

3.6 **Endoscopic retrograde cholangiopancreatography (ERCP):** One of a variety of imaging techniques used to study the biliary tree, whereby an endoscope is passed from the gut via the ampulla of Vater into the biliary ducts.

4. Procedure

4.1. Clinical background:

Bile is a sterile fluid that is made and released by the liver and stored in the gallbladder. Colonization of bile may occur, frequently with a mixture of aerobes and anaerobes originating from the gut. The major cause of biliary infections are Gram negative bacteria (mainly *Escherichia coli*), however Gram positive and anaerobic organisms are also found. Occasionally instrumentation or stenting may lead to infection, which may progress to bacteremia. Fever, previous endoscopic or percutaneous biliary instrumentation, and bilioenteric anastomosis are significant predictors of a positive bile culture. The infection of biliary system can produce significant morbidity and mortality and the prognosis often depends on the presence of biliary

tract obstruction. Biliary infections present as either cholangitis or cholecystitis. The Common organisms isolated from bile are mentioned in the table (1):

Table (1): Common organisms isolated from bile		
Bacteria	*Yeasts	Parasites
Enterobacterales Enterococcus spp. Pseudomonads Bacteroides spp. Clostridium spp. Anaerobes <i>S. aureus</i> Salmonella spp.	<i>Candida albicans</i> , other <i>Candida</i> spp. have been reported.	<i>Ascaris lumbricoides</i> <i>Clonorchis sinensis</i> Opisthorchis spp. <i>Fasciola hepatica</i> <i>Giardia intestinalis</i> Cryptosporidium spp. Microspora

*Rare in normal individuals. They occur in: older patients with malignancy, immunocompromised patients, diabetic patients or in patients receiving antimicrobial treatment for other infections

Note: Other organisms may be isolated and should be given consideration depending on clinical details.

4.2.Pre – analytical stage:

4.2.1. Sample:

- Sample type: Bile fluid (intra-op or from a closed drainage system by aspiration with a needle and syringe) is considered an urgent sample that has to be processed within 2 hrs.
- Sample volume: Minimum requirement is 1 ml.
- Sample stability and storage requirements:
 - Bile samples should be received into appropriate CE marked leakproof containers and place in sealed plastic bags.
 - The optimal time for specimen collection is prior to antimicrobial therapy where possible.
 - Compliance with transport and storage regulations is essential.
 - Specimens should be processed as soon as possible.

- If processing is delayed, refrigeration is preferable to storage at ambient temperature.
- After examination, all specimens are stored for one week (+/-) according to lab storage capacity as additional examinations may be requested during this retention period.

4.2.2. Material:

Table 2: Equipment and supplies

Reagents	Consumables/Supplies	Equipment
Gram stain reagents Agar plates Selenite F Broth Anaerobic jar and anaerobic bag Oil Susceptibility discs	10 µl disposable loops Pasteur pipette	CO ₂ incubator O ₂ incubator (Ambient air) Hot plate Light microscope BSC-II

4.2.3. Safety precaution:

- All specimens need to be treated as potentially infectious.
- Standard procedures for handling of biohazard material must be followed at all times.
- Process bile specimens under BSC-II.
- Processing of diagnostic sample cultures that are assessed to be at higher risk of containing hazard group 3 organisms must be undertaken under Biosafety level-III conditions (BSL-III). Such organisms include Mycobacterium species, Brucella species, Bacillus anthracis, Blastomyces dermatitidis, Histoplasma capsulatum, Coccidioides immitis, etc. ¹

4.2.4. Quality control:

- Check the expiry dates of all media, reagents and stains before use.
- All media, broths, reagents, kits, and stains **MUST** be quality controlled before use and checked for sterility.

- Identification tests should be run with appropriate controls.
- Record the quality control results in the appropriate QC sheet.

4.3. Analytical stage:

4.3.1 Microscopy (Gram stain):

- Label a clean glass-slide with the following details: Bile, lab number, and the date.
- Add one drop of the specimen with a sterile loop to make a thin smear on a clean microscope slide.
- Dry and fix in a hot plate
- Stain with gram stain. Then check under ordinary light microscope.
- Perform an ova/ parasite microscopy ONLY if requested. Using a sterile Pasteur pipette, place one drop of the specimen on a clean glass slide, use with cover slip and scan for parasites.
- Perform TB examination ONLY if requested. ZN stain and TB culture are currently done at CPHL TB Lab-Darsait.

4.3.2 Culture:

- With a sterile Pasteur pipette take 1 drop of the fluid and inoculate it into each agar plate: Blood (Aerobic and Anaerobic) and MacConkey.
- Add few drops of the fluid into selenite F broth.
- Always allow inoculum to dry before spreading to minimise any antibiotic effect which may be present.
- Streak the inoculum using a good streaking technique with a sterile disposable loop. media inoculation should be done in a logical order from least to most selective to avoid the inhibition of organisms by carryover of the selective agent.
 1. Media without inhibitors (Blood)
 2. Indicator media (MacConkey)
 3. Selective media (XLD, Sabouraud (when needed))
- Using forceps, Metronidazole (MTZ) disc is kept between the first and second spread near to the edge (to avoid total inhibition of very susceptible organisms) of anaerobic Blood agar, see figure 1 (clean the forceps by alcohol wipes before and after adding the MTZ disc).




Figure 1: Position of MTZ disk in anaerobic Blood agar plate

- After inoculation, all the plates should be incubated as soon as possible, as per (Table 3):

Table 3: Culture plate selection for Bile specimens						
Clinical/ Gram Stain	Medium	Incubation			Culture read	Significant isolates
		Temp (°C)	Atmosphere	Time (hrs)		
All gram positive and gram negative bacteria	Blood Agar	35-37	5-10 % CO ₂	40- 48	Daily (for 2 days)	Any organism
	Anaerobic Blood Agar	35-37	Anaerobic	48	After 48 h	Anaerobes
	MacConkey	35-37	O ₂	16- 24	Daily (for 2 days)	Enterobacterales
	Selenite F Broth then sub-cultured to XLD	35-37	O ₂	16- 24	Daily (for 2 days)	Salmonella spp. Shigella spp.
Yeast	Sabouraud *	35-37	O ₂	40- 48	Daily (for 2 days)	Candida spp.
* ONLY if yeast seen in Gram stain.						

4.3.3 Identification and Isolation:

- Work-up any growth of potentially pathogenic organisms according to Table (1).
- Minimum level of identification in the laboratory may be set to species level. For fungi, can rely on genus level.

- Any organism considered to be a contaminant/flora may not require identification to species level.
- Work-up a maximum of 3 organisms. Consult in charge technologist or a microbiologist if > 3 pathogens.

4.3.4 Susceptibility Testing:

- As per the national antibiotic susceptibility testing guidelines.

4.4.Post – analytical stage:

4.4.1 Reporting:

4.4.1.1 Appearance: The presence of pus should be noted.

4.4.1.2 Microscopy reporting:

- Quantitate the presence of pus cells (WBC) and organisms (presence of bacteria/yeast / parasite)

Bile Microscopy reporting (LPF)	
Zero cells/organisms	Not seen
<5 / LPF and/or some fields without WBC/ organism	Scanty
5- 10/ LPF	Few
11- 25/ LPF	Moderate
>25 / LPF	Many

4.4.1.3 Culture reporting:

- Preliminary report: 24 hours.
- Final report : 48 - 72 hours
- Report the growth as follows:

Culture result after 48 hours	Reporting comments
Negative report	“No growth after 48 hours of incubation”
Insignificant growth, e.g: CoNS, strept. Viridians..etc)	“No significant growth/ No pathogen isolated”. Add an appropriate comment on possible contamination or overgrowth the specimen is from a collection bag or T-tube)
Any of significant organisms	Report with ID & AST as appropriate (Table 1)

- **Note:**
 - Final positive culture reports are entered in the LIS by the laboratory technician, and then verified and authorized by the medical microbiologist/ senior laboratory technologist.
 - Notify infection control in case of isolation of MDRO's / others as indicated clinically.

5. Responsibilities

5.1. Responsible staff:

- To ensure the adherence to critical result communication procedure
- To facilitate the alternative channels once needed

5.2. Quality manager /officer

- To follow up the implementation of the procedure
- To monitor regularly communication of critical results and raise non-conformance with corrective action once needed.

5.3. All lab staff:

- To adhere to the procedure.
- To document record and release results as recommended
- To report test failures or incident

6. Document History and Version Control

Version	Description	Review Date
1	Initial Release	May 2026

5 References

Title of book/ journal/ art Website	Author	Year of publication	Page
Investigation of Bile.	Public Health England. UK Standards for Microbiology Investigations.	2014	B15 Issue 7.
Inoculation of culture media for bacteriology.	Public Health England. UK Standards for Microbiology Investigations	2017	Q 5 Issue 2.
Biliary tract infections: a guide to drug treatment.	Westphal JF, Brogard JM. Drugs 1999	1999	Drugs 1999;57:81-91. B, II
Clinical factors associated with positive bile cultures during primary percutaneous biliary drainage.	Brody LA, Brown KT, Getrajdman GI, Kannegieter LS, Brown AE, Fong Y et al. JVascIntervRadiol	1998	JVascIntervRadiol 1998;9:5 A, II
Acute cholangitis.	Sinanan MN. InfectDis Clin North Am	1992	InfectDis Clin North 1992;6:571-99. B, II
Association of preoperative biliary stenting with increased postoperative infectious complications in proximal cholangiocarcinoma.	Hochwald SN, Burke EC, Jarn WR, Fong Y, Blumgart ArchSurg	1999	ArchSurg 1999;134:261-6. B,
Management of ascending cholangitis.	Kinney TP. GastrointestEndoscClinNAM	2007	GastrointestEndoscClinNAM 2007;17:289-306, vi. B, II
Clinical determinants for the recovery of fungal and mezlocillin-resistant pathogens from	Ehrenstein BP, Salamon L, Linde HJ, Messmann H, Scholmerich J, Gluck T. Clinical Infectious Diseases	2002	Clinical Infectious Dis 2002;34:902-8. B, II

specimens.			
Fungal cholecystitis and AIDS	Cervia JS, Murray HW. J InfecDis	1990	J InfectDis 1990;161:358. A, I