





Ministry of Health

| | |
|--------------------------------|---|
| Document Title | Stool culture SOP |
| Document Type | Procedure |
| Directorate/Institution | Directorate General of Specialized Medical Care |
| Targeted Group | Medical laboratories |
| Document Author | Dr. Noora Al Busaidi |
| Designation | Medical officer |
| Document Reviewer | Microbiology documents development team |
| Designation | Microbiology documents development team |
| Release Date | May 2023 |
| Review Frequency | Three Years |

| Validated by | | Approved by | |
|---------------------|---|--------------------|---|
| Name | Dr. Muna Habib | Name | Dr. Badryah Al Rashidi |
| Designation | Director Department Development & Control (DGQAC) | Designation | Director General of Primary Health Care |
| Signature |  | Signature |  |
| Date | May 2023 | Date | June 2023 |

Contents Table:

| | |
|--|----|
| Acronyms:..... | 4 |
| .1 Purpose..... | 5 |
| 2. Scope..... | 5 |
| 3. Definitions..... | 5 |
| 4. Procedure..... | 5 |
| 5. Responsibilities..... | 9 |
| 6. Document History and Version Control..... | 11 |
| 7. References:..... | 11 |
| 8. Annexes..... | 12 |

Acknowledgment

The diagnostic laboratories services at the Directorate General of Specialized Medical Care (DGSMC) at Ministry of Health (MOH) would like to thank and appreciate the great effort of the Microbiology documents development team. Participated and contributed personnel are:

| Member name | Institution | Designation |
|--------------------------------|--|---|
| Dr.Mahmoud Al Subhi | Rustaque Hospital | Team Leader Consultant microbiologist |
| Ms. Zainab Al Hadhrami | Directorate General of Specialized Medical Care | Team Coordinator Head of Quality unit in Diagnostic laboratories services |
| Ms. Saleh Al Shukairi | Ibra Hospital | Senior technologist specialist A |
| Dr. Hanaa Al Auraimi | Royal Police of Oman Hospital | Consultant microbiologist |
| Dr. Nawal AL Kindi | Khoula Hospital | Consultant microbiologist |
| Dr. Al Warith Al Kharousi | Nizwa Hospital | Consultant microbiologist |
| Dr. Abdulrahman Al Mahrouqi | Ibri Hospital | Specialist microbiologist pathologist |
| Dr. Nada Al Tamtami | Al Massara Hospital | Consultant microbiologist |
| Dr. Wafaa Al Tamtami | Armed Forces Hospital | Consultant microbiologist |

Acronyms:

| | |
|------|---|
| BA | Blood agar |
| MAC | Mackonky plate |
| SMAC | Sorbitol mackonky plate |
| TCBS | Thiosulphate Citrate Bile Salt Sucrose Agar |
| XLD | Xylose Lysine Deoxycholate agar |
| APB | Alkaline Phosphate Broth |
| H&S | Health and Safety |
| ID | Identification |
| IQC | Internal Quality Control |
| SOP | Standard operating procedure |
| NICU | Neonatal intensive care unit |
| KIA | Kliger's Iron Agar |

1. Purpose

This document describes the procedure for stool sample culture.

2. Scope

This document is applicable for all medical laboratories under MOH and other collaborative governmental and non-governmental health institutions.

3. Definitions

3.1 Diarrhea: defined as stools of a loose or liquid consistency, occurring more frequently than is normal for the individual and usually this is at least three or more instances in a 24-hour period.

4. Procedure

4.1. Clinical background:

Acute infectious diarrhea may be caused by bacteria, viruses or protozoa. The laboratory will routinely look for the most common causes of bacterial diarrhea. This includes *Campylobacter*, *Salmonella*, *Shigella*, *E.coli* O157:H7, and *Yersinia*. Less common causes such as *Aeromonas*, *Plesiomonas*, and *Vibrio* will be cultured only when clinically indicated or specifically requested.

There are three types of diarrhea:

1. Acute watery diarrhea: This is defined as diarrhoea not exceeding 14 days of duration and typically with frequency of 3 or more episodes a day. Organisms implicated include norovirus, rotavirus, adenovirus, *Salmonella* species, *Campylobacter* species and *C. difficile*. Other organisms that may be considered include the *Vibrio* species including *Vibrio cholerae* (the causative agent of cholera).
2. Acute bloody diarrhoea: This is a sudden onset of diarrhoea where frank blood is present. Dysentery is an acute infectious gastroenteritis characterized by loose stools with blood and mucus. Organisms implicated in acute bloody diarrhoea include *Campylobacter* species, *Shigella* species, Shiga toxin producing *Escherichia coli* (STEC) including serogroup O157, *Salmonella* species, *Entamoeba histolytica*.
3. Persistent and chronic diarrhoea: persistent diarrhoea is defined as diarrhoea of >14 days but fewer than 30 days in duration. It should be noted that viruses (eg norovirus) and bacteria (*Salmonella*, *Shigella* and *Campylobacter* species) can be the cause of

persistent diarrhoea in patients who are immunocompromised. Diarrhoea that lasts >30 days is referred to as “chronic”. Organisms implicated are predominantly parasites – *Giardia* species, *Cryptosporidium* species, *Cyclospora cayetanensis* and *Microsporidia* species. *Chlamydia trachomatis* infection including lymphogranuloma venereum (LGV) and *Neisseria gonorrhoeae* infection can present atypically as sexually transmitted infectious colitis in individuals who have unprotected anorectal sex.

4.2.Pre – analytical stage:

4.2.1. Sample:

- A single stool specimen for routine culture for enteric pathogens should be collected into a leak-proof container and transported to the laboratory as soon as possible after collection, Specimens should be processed upon arrival to the laboratory immediately, because the number of important pathogens such as *Shigella* species may not survive the pH changes that occur in feces specimens which are not promptly delivered to the laboratory even if refrigerated.
- Samples from hospitalized patients > 72 hours are not accepted (three day rule) and reported as hospital acquired diarrhea and released with comment: Diarrheal sample specimens from patients who have been admitted for more than 3 days are not routinely processed for stool culture, if clinically indicated please call microbiologist.

4.2.2. Material:

| Reagents | Supplies | Equipment |
|--|---|---|
| Selenite F MAC, BA, XLD, Sorbitol MAC, CAMPY, TCBS 70% alcohol reagent. KIA TSI slants Distilled water. 6.5% Normal saline. Campy gas kit generator salmonella, sheigella, E.coli antisera | Disposable loops. wooden sticks, Racks. Plastic streaking loops. Universal plastic containers. Waste disposable, Sterility disposable bags. | 43°C incubator O ₂ incubator Fume hood Jars |

- All the antisera has to be kept in fridge (2-8 C).

4.2.3. Safety precaution:

- All specimens need to be treated as potentially infectious. Standard procedures for handling of biohazard material must be followed at all times. Universal Precautions must be practiced at all stages of these procedures.
- Specimen processing in general requires containment level 2 for routine work. Containment level 3 if the following organisms are suspected from clinical information or laboratory findings: *Salmonella Typhi*, *Salmonella Paratyphi A, B, C*, *E.Coli 0157:H7*, *Shigella dysenteriae*, *Toxigenic V.cholerae*.

4.2.4. Quality control:

- Check the expiry dates of all media, reagents and stains before use.

- All media, reagents, kits, and stains MUST be quality controlled before use.
- Identification tests should be run with appropriate controls.
- Record the quality control results in the appropriate QC sheet.
- Monitor the 43 ° incubator temperature daily.

4.3. Analytical stage:

4.3.1 Direct exam: refer to stool microscopy SOP.

4.3.2 Culture set up:

1. Take the received fecal sample to the fume hood.
2. Prepare the appropriate plates according to the following table:

| Media | Incubation |
|--------------------------------|---------------------------------------|
| MacConkey Agar (MAC) | O ₂ , 35°C x 18 – 24 hours |
| XLD | O ₂ , 35°C x 18 - 24 hours |
| MacConkey Sorbitol Agar (SMAC) | O ₂ , 35°C x 18 - 24 hours |
| Campylobacter Agar (CAMPY) | Campy Jar, 42°C x 48 hours |
| Selenite F Broth (SEL) | O ₂ , 35°C x 12-18 hours |

- If *Yersinia* is requested or patient is >1 month–12 years old (except for NICU) OR If history of mesenteric adenitis, terminal ileitis add:
Yersinia agar (if available) O₂, 30°C x 48 hours OR incubate MAC for 48hrs at room temperature
- If *Vibrio* is requested, or “if there is a history of shellfish ingestion or if there is a history of travel outside to areas where cholera/Vibrio associated diarrhoea is endemic”, add: TCBS O₂, 35°C x 18 - 24 hours
APB O₂, 35°C x 5 - 8 hours
- If *Plesiomonas* or *Aeromonas* is requested , or “if there is a history of sea food ingestion or if there is a history of travel outside to areas add: Blood Agar (BA) O₂, 35°C x 24 hours

4.3.3 Isolation and Identification:

1. After 24 hrs of incubation: Identify any growth of *Salmonella*, *Shigella*, *E. coli* O157, or *Yersinia* by examining the MAC, XLD, and SMAC at 24 hr.
2. After overnight incubation, subculture Selenite broth onto XLD which is incubated in O₂ at 35°C x 24 hr. Workup as indicated above.
3. Examine MAC again at 48hr for (if *Yersinia* was requested or indicated by age group (>1 month–12 years old)
4. Examine Campy plates at 48 hr and identify any growth of *Campylobacter*.
5. See Identification schemes /flow charts for details. These, however, are meant as guides only.
6. Identify other pathogens only if specifically requested.

4.3.4 Screening for Salmonella and Shigella:

1. Examine MAC and DCA/XLD for suspected colonies:
2. Colourless or transparent NLF or very pale pink like LLF(*Shigella*) on MAC and oxidase negative .
 - Pink to red with or without black centres on XLD.
 - Colourless with or without black centres on DCA.
3. Pick one colony of each suspect morphotype and inoculate a KIA, peptone water. Incubate each according to protocol at 35°C in O₂ and add BA as purity plate.
4. Determine if the organism is a potential enteric pathogen according to results of Table 1.
5. Suggestive screening results require serologic confirmation and SET UP for available automated or manual identification system.
6. Set up of ID and serologic confirmation can be done directly from growth on the blood agar plate.
7. If biochemical tests indicate salmonella, send the isolate directly for susceptibility.

Table 1. Characteristic reactions of potential enteric pathogens and follow-up

| KIA/gas/H ₂ S | Urea | Indole | Motility | Follow-up |
|--|------|-------------------|--|---|
| K/A, + gas, + H ₂ S (degree of gas and H ₂ S may vary) ^{1,2} | - | - | + | Salmonella Anti-sera , API20E, AST and Send to PHL |
| K/A, no gas H ₂ S neg. ³ | - | - / + | - | Shigella Anti-sera , API20E, AST and Send to PHL |
| A/A, +/- gas | +/- | + | +/- | No further testing |
| K/K | +/- | Not applicable | +/- | No further testing |
| K/A or A/A No gas, no H ₂ S or “rainbow” | + | + | Positive at room temp. Negative at 35°C | Pinpoint on MAC Yersinia API20E, AST and Send to PHL |

Note: ¹*S. Typhi* may only produce small amounts of gas and/or H₂S.

²*S. Paratyphi A* may occasionally produce H₂S weakly.

³*Shigella Flexner* (type 6) may produce a small amount of gas.

Table 2. Characteristic reactions of potential enteric pathogens

| Bacteria | Urea | Slant | Butt | H₂S | Gas | Motility | Indole | Oxidase |
|--------------------------|-------------|--------------|-------------|-----------------------|------------|-----------------|---------------|----------------|
| E. coli | - | A | A | - | + | (+) | d | - |
| E. coli - Alk. Despair | - | K | A | - | - | - | + | - |
| Klebsiella | + | A | A | - | + | - | d | - |
| Enterobacter | - | A | A | - | + | + | - | - |
| Citrobacter | dw | d | A | d | + | + | - | - |
| Salmonella | - | K | A | dw | d | + | - | - |
| S. typhi | - | K | A | + | - | + | - | - |
| S. paratyphi A | - | K | A | - | + | + | - | - |
| S. arizonae | - | d | A | + | + | + | - | - |
| Shigelladysenteriae | - | K | A | - | - | - | d | - |
| S. flexneri | - | K | A | - | -* | - | d | - |
| S. boydii | - | K | A | - | -** | - | d | - |
| S. sonnei | - | K | A | - | - | - | - | - |
| Proteus vulgaris | + | K | A | + | d | + | + | - |
| P. marabilis | + | K | A | + | + | + | - | - |
| P. morganii | + | K | A | - | d | d | - | - |
| P. rettgeri | + | K | A | - | d | + | + | - |
| Providenciaalcalifaciens | - | K | A | - | d | + | + | - |

| | | | | | | | | |
|-------------------------|---|---|-----|---|---|---|---|---|
| P. stuartii | - | K | A | - | - | + | + | - |
| Yersinia enterocolitica | + | K | A | - | - | V | d | |
| Aeromonas | - | K | A | - | + | + | + | + |
| Plesiomonas | - | K | A | - | - | + | + | + |
| Vibrio | - | K | A | - | - | + | + | + |
| Alkaligenes | - | K | N/K | - | - | + | - | + |
| Edwardsiella | - | K | A | + | + | + | + | - |
| Hafnia | - | d | A | - | + | + | - | - |

Symbols:

K = alkaline (red) reaction

A = acid (yellow) reaction

N = neutral reaction

+ = positive reaction

w = weak or delayed reaction

V = variable reaction (Y. enterocolitica is motile at 25° C but non-motile at 37°C.

d = 10-85% positive

* Some S. flexneri serotype 6 gas (+)

** Serotype 13 and 14 gas (+).

4.3.5 Susceptibility testing:

After primary identification, follow the national antimicrobial susceptibility guidelines for all isolated enteropathogens.

4.4.Post – analytical stage:

4.4.1. Reporting :

- Negative Report: release the negative (non-significant) stool culture results as: *Salmonella*, *Shigella*, *Campylobacter* are NOT isolated. If specially requested for *Yersinia* or *E.coli O157*, include them pathogens to the comment.
- Positive Report: release the “Name of pathogen” ISOLATED.
- There is no need for quantification as any amount is significant.

- There is no need to comment on the absence of the other pathogens normally screened for.
- In case of any campy QC failure, release the stool culture report with the comment “campylobacter was not tested, if patient is suspected to have campylobacter infection please submit another fresh stool sample”.
- Follow the procedure for reporting of communicable enteric pathogens.
- Notify microbiologist/ senior of all *Shigella*, *E. coli O157*, *Salmonella Typhi* or *Paratyphi*, or *Vibrio cholera*.
- Refer to critical results notification S.O.P.

5. Responsibilities

5.1.Responsible staff:

- To ensure the adherence to critical result communication procedure
- To facilitate the alternative channels once needed

5.2.Quality manager /officer

- To follow up the implementation of the procedure
- To monitor regularly communication of critical results and raise non-conformance with corrective action once needed.

5.3.All lab staff:

- To adhere to the procedure.
- To document record and release results as recommended
- To report test failures or incident

6. Document History and Version Control


| Version | Description | Review Date |
|---------|-----------------|-------------|
| 1 | Initial Release | May 2026 |
| | | |
| | | |
| | | |

7. References


| Title of book/ journal/ articles/ Website | Author | Year of publication | Page |
|--|--------|---------------------|---------------|
| 1. Church et al. Am J ClinPathol Siegel 103 | | 1995; | 149- 53. |
| 2. Al. JAMA; 263 | | 1990 | 979- 82. |
| 3. Brazier JS. JAC; 41, Suppl. | | 1998 | C: 29- 40. |
| 4. UK SMI S 7 | | 2020 | |

8. Annexes


A. Stool culture algorithm




B. Screening for salmonella and shigella




C. Identification of salmonella




C. identification of salmonella continued




D Identification of shigella




E. Identification of campylobacter



F. Identification of E.coli O157:h7



G. Screening primary MacConkey Agar for *Yersinia*



H. Screening primary BA plate for *Vibrio*, *Aeromonas* & *Plesiomonas*

