

# Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs  
and Drug Control  
MUSCAT



سِلاطِنَا اَمَانَا  
وَزَارَةُ الصِّحَّةِ  
المديرية العامة للصحة  
والرقابة الدوائية  
مسقط

To:

**THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES**

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

**ALL PRIVATE PHARMACIES & DRUG STORES**

After Compliments,

Please find attached our Circular No...<sup>79/2020</sup> dated ...<sup>6/14/2020</sup> Regarding Recall of Low pressure tube, CO2 from Karl Storz GmbH

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information

# Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs  
and Drug Control  
MUSCAT



سلطنة عمان  
وزارة الصحة  
المديرية العامة للصيدلانية  
والرقابة الدوائية  
مسقط

Circular No. 79 / 2020

12-08-1441 H

6-04-2020

Ref: 48/2020

## Recall of Low pressure tube, CO2 from Karl Storz GmbH

Source of Recall	NCMDR National Centre for Medical Devices Reporting, SFDA <a href="https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=15100">https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=15100</a>
Product	Low pressure tube, CO2- Physical therapy ,insufflation technique
Manufacturer	Karl Storz GmbH & Co KG
Local Agent	Mustafa Sultan Science & Industry
The affected products	Article numbers: UI400, UI500, 26432020-1, 26430520-1, UI001; UI002; UI003; Serial numbers: all
Reason for Recall	Karl Storz has been informed that the low pressure CO2 tubes stated above may contain an incorrect component in the connecting plug. The diameter of the component deviates from the specified dimension
Action	1. Kindly check your stock, contact your local agent for remedial action 2. Locate all potentially affected products and check their inner diamete 3. Effected products must be returned to Karl Storz Contact your local agent to arrange for a replacement
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical Device to Director of Medical Device Control contact E-mail <a href="mailto:Med-device@moh.gov.om">Med-device@moh.gov.om</a>

Directorate General of Pharmaceutical Affairs & Drug Control  
Sultanate of Oman

Dr. Mohammed Hamdan Al Rubaie

DIRECTOR GENERAL

