



To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES
Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)
Director General of Engineering Affairs, MOH
Director General of Royal Hospital
Director General of Khoula Hospital
Director General of Medical Supplies (MOH)
Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)
Hospital Director (Al Nahda Hospital)
Hospital Director (Al Massara Hospital)
The Head of Medical Services in SQU Hospital
The Head of Medical Services in Royal Oman Police
The Head of Medical Services in Ministry of Defence
The Head of Medical Services in The Diwan
The Head of Medical Services in The Sultan's Special Force
The Head of Medical Services in Internal Security Services
The Head of Medical Services in Petroleum Development of Oman
The Head of Medical Services in LNG Oman
ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No. 99..... dated 22/5/22 Regarding GHC recall of Sevoflurane vaporizers for Flow Family Anesthesia systems from (mrf: GETINGE).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information



Circular No. 99 / 2022


بنقدم
Forward
with Confidence

رؤية عمان
2040
Oman

23 -10-1443 H

22-05-2022

Recall of Sevoflurane vaporizers for Flow Family Anesthesia systems from GETINGE.

Source	GHC- Gulf Health Council.						
Product	Sevoflurane vaporizers for Flow Family Anesthesia systems.						
Manufacturer	GETINGE.						
The Affected Products	<table border="1"><thead><tr><th>Affected vaporizer</th><th>Replacement vaporizer</th></tr></thead><tbody><tr><td>6886601, Vaporizer Sevoflurane, Maquet filling</td><td>6682282 Vaporizer Sevoflurane, Maquet filling</td></tr><tr><td>6887523, Vaporizer Sevoflurane, SAFE-T-SEAL</td><td>6887135, Vaporizer Sevoflurane, SAFE-T-SEAL</td></tr></tbody></table>	Affected vaporizer	Replacement vaporizer	6886601, Vaporizer Sevoflurane, Maquet filling	6682282 Vaporizer Sevoflurane, Maquet filling	6887523, Vaporizer Sevoflurane, SAFE-T-SEAL	6887135, Vaporizer Sevoflurane, SAFE-T-SEAL
Affected vaporizer	Replacement vaporizer						
6886601, Vaporizer Sevoflurane, Maquet filling	6682282 Vaporizer Sevoflurane, Maquet filling						
6887523, Vaporizer Sevoflurane, SAFE-T-SEAL	6887135, Vaporizer Sevoflurane, SAFE-T-SEAL						
Local Agent	Waleed Pharmacy & Stores LLC.						
Reason	Hydrogen fluoride has been confirmed by chemical analysis. Potential harms may include irritation of respiratory tract and in worst case may lead to lung edema and/or severe hypocalcemia.						
Action	1. All HCP's should Immediately stop using the affected products and contact the supplier for replacement. 2. Contact the local agent for remedial action.						
Product Picture							
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control contact E-mail: Med-device@moh.gov.om						

Dr. Mohammed Hamdan Al-Rubaie
DIRECTOR GENERAL

