


Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs
and Drug Control

MUSCAT



سِلاطِنَا مَومَانَا
وَزَارَةَ الصِّحَّةِ
وَالرِّقَابَةِ وَالدُّوَالِيَّةِ
مَسْقَط

To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No. 117 dated 17/06/20 Regarding Recall of safeCLINITUBES Blood Samplers (Mfr: Radiometer)

Copy to:


- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information

Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs
and Drug Control

MUSCAT




سلطنة عمان
وزارة الصحة
المديرية العامة للصحة
والدواء والرقابة الدوائية
مسقط

Circular No. 117 / 2020

25 -Shawal-1441 H

17 -June-2020

Recall of safeCLINITUBES Blood Samplers From Radiometer.

| | |
|-----------------------|---|
| Source | NCMDR (National Center for Medical Devices Reporting). https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&rid=15121 |
| Product | safeCLINITUBES Blood Samplers. Tubes, Blood Collection, Capillary/Microsample. |
| Manufacturer | Radiometer. |
| The affected products | Order number : (Lot number range) - 942-892: (R2425 - R2524) - 942-898: (R1431 - R1495) - 942-969: (R0291 - R0301) |
| Reason | May Have Insufficient Aspiration Capabilities, Potentially Resulting in Lost Samples. |
| Action | 1. Cease using the safeCLINITUBES Blood Samplers from the affected lots. 2. Check your inventory of safeCLINITUBES Blood Samplers from the affected lots. 3. Contact the local agent for replacement of affected product if available. |
| Product image |  |
| comments | Healthcare professionals are encouraged to report any adverse events suspected to be associated with the above device or any other medical device to Department of Medical Device Control contact E-mail Med-device@moh.gov.om |

Dr. Mohammed Hamdan Al Rubaie

DIRECTOR GENERAL

Directorate General of Pharmaceutical Affairs & Drug Control
Sultanate of Oman

