


Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs
and Drug Control

MUSCAT



سِيَّاطِنَا جَمَانَا
وَزَارَةُ الصِّحَّةِ
المديرية العامة للصحة
والرقابة الدوائية
مسقط

To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No. 152 dated 01/09/2012 Regarding NCMDR Recall of Liaison® Measles from (mfr: DiaSorin Inc).

Copy to:


- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information

Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs
and Drug Control

MUSCAT




سلطنة عمان
وزارة الصحة
المديرية العامة للأدوية
والرقابة الدوائية
مسقط

Circular No. 162/2020

12 -1-1442 H

01 -09-2020

Recall of Liaison® Measles from DiaSorin Inc.

Source	NCMDR- National Centre Medical Device Reporting https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&rid=15298
Product	Liaison® Measles.
Description	IVD .
Manufacturer	DiaSorin Inc.
The affected products	Catalog Number: 318820 Lot Number: Lot 174026 - Expiry date: 2021/08/08
Reason	DiaSorin S.p.a. internal testing confirmed, for the lot in object, a reduction of index value of samples. Such decrease could lead to the risk of False Negative results classification in a limited number of tests on samples having low Measles virus IgM concentration, closer to the assay Cut-Off (i.e. Index 1.0).
Action	1. Please, stop using this lot of Liaison® Measles IgM kit on Liaison Analyzer Family and destroy any remaining inventory. 2. Contact the local agent for remedial action.
Product image	
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control contact E-mail: Med-device@moh.gov.om

Dr. Mohammed Hamdan Al Rubaie

DIRECTOR GENERAL

