

To:

**THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES**  
**Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)**  
**Director General of Engineering Affairs, MOH**  
**Director General of Royal Hospital**  
**Director General of Khoula Hospital**  
**Director General of Medical Supplies (MOH)**  
**Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)**  
**Hospital Director (Al Nahda Hospital)**  
**Hospital Director (Al Massara Hospital)**  
**The Head of Medical Services in SQU Hospital**  
**The Head of Medical Services in Royal Oman Police**  
**The Head of Medical Services in Ministry of Defence**  
**The Head of Medical Services in The Diwan**  
**The Head of Medical Services in The Sultan's Special Force**  
**The Head of Medical Services in Internal Security Services**  
**The Head of Medical Services in Petroleum Development of Oman**  
**The Head of Medical Services in LNG Oman**  
**ALL PRIVATE PHARMACIES & DRUG STORES**


After Compliments,

Please find attached our Circular No. **1.62** dated **13/9/22** Regarding NCMDR Recall of Rainbow Super DCI-mini sensor (SpCO, SpHb, SpMet, and SpO2) from (mfr: Masimo Corporation.).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information





Circular No. 168/ 2022


تقدم بثقة  
Moving Forward  
with Confidence



17 -02-1444 H

13 -09-2022

**Recall of Rainbow Super DCI-mini sensor (SpCO, SpHb, SpMet, and SpO2) from Masimo Corporation.**

Source	NCMDR- National Center for Medical Devices Reporting- SFDA <a href="https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=17256">https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=4&amp;rid=17256</a>
Product	Rainbow Super DCI-mini sensor (SpCO, SpHb, SpMet, and SpO2)
Description	Spot check sensor.
Manufacturer	Masimo Corporation.
Local agent	Mustafa Sultan Science & Industry Co.L.L.C.
The affected products	Part # 4309 LOT 20JCP and Part # 4311 LOT 20FFP.
Reason	May provide higher-than expected SpCO readings on patients with normal COHb levels.
Action	1. Return of the affected part and lot number(s) for replacement. 2. Contact the local agent for remedial action.
Product Image	
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: <a href="mailto:Med-device@moh.gov.om">Med-device@moh.gov.om</a>

Dr. Mohammed Hamdan Al Rubaie  
DIRECTOR GENERAL

