



بمقدم بثقة
Moving Forward
with Confidence



To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No. 76 dated 2014/22 Regarding NCMDR FSN of Compella Bariatric Bed Systems and Therapy Air Supply Units from (mrf: Hill Rom Inc).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information



PADC
المديرية العامة للصيدلة والرقابة الدوائية
Directorate General of Pharmaceutical
Affairs & Drug Control



ص.ب: 393 مسقط - الرمز البريدي: 100 - هاتف: 22357111 - فاكس: 22358489

P.O. Box: 393 Muscat - Postal Code : 100 - Tel: 22357111 - Fax: 22358489

dgpa_dc Email: dg-padc@moh.gov.om





Circular No. 76 / 2022

18-09-1443 H

20-04-2022

Field Safety Notice of Compella Bariatric Bed Systems and Therapy Air Supply Units from Hill Rom Inc

Source	NCMDR-National Center for Medical Device Reporting https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=6&rid=16093
Product	Compella Bariatric Bed Systems and Therapy Air Supply Units.
Description	Beds, Electric, Bariatric.
Manufacturer	Hill Rom Inc.
Local Agent	Taiba Medserv LLC.
The affected products	Bariatric Bed Systems Model No. P7800A Manufacturing Dates: All beds manufactured. Therapy Air Supply Unit Model No. P7810A Manufacturing Dates: All Air Supply units manufactured
Reason	Powers cord can become damaged because of improper storage that exposed AC power conductors due to cut or damaged power cords.
Action	1. Please follow Hill Rom recommendations to mitigate the risk and please share with potential users in your organization as a reminder to follow the instructions for proper handling techniques. 2. If a power cord is damaged, will be replaced. 3. Contact the local agent for remedial action.
Product Picture	  <p>Fig. 1- Compella Cord Hooks. Fig. 2- Power Cord Magnets</p>
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control contact E-mail: Med-device@moh.gov.om

Dr. Mohammed Hamdan Al Rubaie
DIRECTOR GENERAL

