Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs and Drug Control

MUSCAT





To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES

Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)

Director General of Engineering Affairs, MOH

Director General of Royal Hospital

Director General of Khoula Hospital

Director General of Medical Supplies (MOH)

Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)

Hospital Director (Al Nahda Hospital)

Hospital Director (Al Massara Hospital)

The Head of Medical Services in SQU Hospital

The Head of Medical Services in Royal Oman Police

The Head of Medical Services in Ministry of Defence

The Head of Medical Services in The Diwan

The Head of Medical Services in The Sultan's Special Force

The Head of Medical Services in Internal Security Services

The Head of Medical Services in Petroleum Development of Oman

The Head of Medical Services in LNG Oman

ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information

Sultanate of Oman

Ministry of Health

Directorate General of Pharmaceutical Affairs and Drug Control

MUSCAT



ميت المكنى بين الى في الى في المنظمة المنظمة

Circular No. of /2021

03 -06-1442 H

1 -01-2021

Recall of Hasson Trocar from Unimax Medical Systems Inc.

Source	NCMDR- National Centre for Medical Devices Reporting- SFDA https://ncmdr.sfda.gov.sa/Secure/CA/CaViewRecall.aspx?caid=6&rid=15505
Product	Hasson Trocar.
Description	Medical instruments for use in humans - visceral surgery.
Manufacturer	Unimax Medical Systems Inc.
Local agent	Medical & Scientific Supplies LLC.
The affected products	Hasson trocar 12mm Transp. Cannula 10cm, Component code: 899307, 07323190273140, 6292002020
Reason	The protective flanges come away from trocar cannula.
Action	 Identify and quarantine affected devices. Contact the local agent for remedial action.
Product image	
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through an E-mail: Med-device@moh.gov.om

Dr. Mohammed Hamdan Al Rubaie

Director General



XX.