



To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES  
Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)  
Director General of Engineering Affairs, MOH  
Director General of Royal Hospital  
Director General of Khoula Hospital  
Director General of Medical Supplies (MOH)  
Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)  
Hospital Director (Al Nahda Hospital)  
Hospital Director (Al Massara Hospital)  
The Head of Medical Services in SQU Hospital  
The Head of Medical Services in Royal Oman Police  
The Head of Medical Services in Ministry of Defence  
The Head of Medical Services in The Diwan  
The Head of Medical Services in The Sultan's Special Force  
The Head of Medical Services in Internal Security Services  
The Head of Medical Services in Petroleum Development of Oman  
The Head of Medical Services in LNG Oman  
ALL PRIVATE PHARMACIES & DRUG STORES


After Compliments,

Please find attached our Circular No. 171 dated 18/9/2022 Regarding NCMDR Recall of IH-QC8 from (mfr: DiaMed GmbH).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information





**Circular No. 171 / 2022**


بنفسه بثقة  
Moving Forward  
with Confidence



21 -02-1444 H

18 -09-2022

**Recall of IH-QC8 from DiaMed GmbH.**

Source	NCMDR- National Center for Medical Devices Reporting- SFDA <a href="https://ncmdr.sfda.gov.sa/Secure/CA/CAViewRecall.aspx?caid=4&amp;rid=17271">https://ncmdr.sfda.gov.sa/Secure/CA/CAViewRecall.aspx?caid=4&amp;rid=17271</a>
Product	IH-QC8.
Description	IVD.
Manufacturer	DiaMed GmbH.
Local agent	Bahwan Healthcare Center LLC.
The affected products	Catalog No. : 009328 IHD Lot n°: 08780 82 1 SAP Lot n°: 746258821 UDI: (01)07611969231560, (17)220919(10)746258821 Expiry Date: 022-09-19
Reason	Tubes from the impacted lot may present clots of red blood cells and turbid supernatant.
Action	1. Stop using the impacted lot. 2. Use another lot of IH-QC8. (Lot 08780 84 1 is now available) 3. Contact the local agent for remedial action.
Product image	
comments	Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: <a href="mailto:Med-device@moh.gov.om">Med-device@moh.gov.om</a>

**Dr. Mohammed Hamdan Al Rubaie**

**Director General**

