



بتقدّر بثقة
Moving Forward
with Confidence

رؤية عمان
2040
Vision

To:

THE DIRECTOR GENERAL OF HEALTH SERVICES IN ALL GOVERNORATES
Commanding Officer, Armed Forces Hospital (Al Khoudh & Salalah)
Director General of Engineering Affairs, MOH
Director General of Royal Hospital
Director General of Khoula Hospital
Director General of Medical Supplies (MOH)
Director General of Pvt. Health Est. Affairs (to kindly arrange distribution to all Pvt. Hospitals)
Hospital Director (Al Nahda Hospital)
Hospital Director (Al Massara Hospital)
The Head of Medical Services in SQU Hospital
The Head of Medical Services in Royal Oman Police
The Head of Medical Services in Ministry of Defence
The Head of Medical Services in The Diwan
The Head of Medical Services in The Sultan's Special Force
The Head of Medical Services in Internal Security Services
The Head of Medical Services in Petroleum Development of Oman
The Head of Medical Services in LNG Oman
ALL PRIVATE PHARMACIES & DRUG STORES

After Compliments,

Please find attached our Circular No 108 dated 29/5/23 Regarding NCMDR Recall of of 3M™ Attest™ Steam Chemical Integrators from (mfr: 3M Health Care).

Copy to:

- Director, Office of H.E. The Undersecretary for Health Affairs
- Director of Medical Device Control, DGPA&DC
- Director of Pharmacovigilance & Drug Information Dept, DGPA&DC
- Director of Drug Control Department, DGPA&DC
- Director of Pharmaceutical Licensing Department, DGPA&DC
- Director of Central Quality Control Lab., DGPA&DC
- Supdt. of Central Drug Information



PADC
المديرية العامة للصيدلة والرقابة الدوائية
Directorate General of Pharmaceutical
Affairs & Drug Control



ص.ب: ٣٩٣ مسقط - الرمز البريدي: ١٠٠ - هاتف: ٢٢٣٥٧١١١ - فاكس: ٢٢٣٥٤٤٨٩

P.O. Box: 393 Muscat - Postal Code: 100 - Tel: 22357111 - Fax: 22358489

dgpa_dc Email: dg-padc@moh.gov.om



Circular No. 108 / 2023

09 -11-1444 H

29 -05-2023

بنقدم
Moving Forward
with Confidence

رؤية عمان
2040
Oman

Recall of 3M™ Attest™ Steam Chemical Integrators from 3M Health Care.

| | |
|-----------------------|--|
| Source | NCMDR - National Center Medical Device Reporting- SFDA. https://ncmdr.sfda.gov.sa/Secure/CA/CAViewRecall.aspx?caid=4&rid=19525 |
| Product | 3M™ Attest™ Steam Chemical Integrators. |
| Description | Sterilization Monitor. |
| Manufacturer | 3M Health Care. |
| Local agent | Waleed Pharmacy & Stores. |
| The affected products | Catalog Number: 1243A Lot Code: FD122025 Date of Manufacture: 27 Dec 2022 Expiration Date: 26 Dec 2025 |
| Reason | Increased potential to leak ink during a sterilization cycle and may provide an unclear result of the acceptability of the sterilization cycle. |
| Action | 1. Please identify the affected product listed above and immediately cease use of affected lot. 2. Please discard all remaining affected product. 3. Contact the local agent for remedial action. |
| comments | Healthcare professionals are encouraged to report any adverse events Suspected to be associated with the above device or any other medical device to Department of Medical Device Control through the E-mail: Med-device@moh.gov.om |

Dr. Mohammed Hamdan Al Rubaie
Director General

